



MEDICAL · SPECIALTY MARKETS · CONSTRUCTION


Embroidery Options

- * Custom embroidery and stock designs available.
- * A logo or design specific to your facility will further enhance the product and clearly identify it.
- * For more options, please refer to our main Medical Catalogue.



Thyroid Shields

- * .30mm (for dental application only) or .50mm L.E.
- * Available in all fabric colors.
- * Can be ordered Loose or Integrated (sewn on).
- * For integrated shields, please see apron styles.



CODE (.30mm)	CODE (.50mm)	DESCRIPTION
TSRL-3	TSRL-5	Standard Thyroid Shield Loose
TSML-3	TSML-5	Mammography Style Thyroid Shield Loose
TSDL-3	TSDL-5	Dental Thyroid Shield Loose

Technician Apron

- * Full length: Protects from chest to knees.
- * 1" straps cross the back and tie securely at the front.
- * Available in XS-XXL, with or without integrated (sewn-on) thyroid collar.
- * Easy on and off convenience.
- * Optional no-charge chest pocket.
- * Available in nylon or vinyl fabric options.

CODE	APRON DESCRIPTION
TA-3	Technician apron tie strap .30mm LE
TA-5	Technician apron tie strap .50mm LE



Tie-strap back

	Shoulder Width	Max. Chest/Hip	Apron Length
X-Small	14"	34"	34"
Small	14.5"	38"	35"
Medium	16"	42"	37"
Large	18.5"	44"	39"
X-Large	20.5"	46"	41"
XX-Large	21.5"	50"	43"



MDEL 2312





Patient Apron (Adult/Child)

- * Available with or without thyroid collar.
- * Hanging loops allow for proper storage.
- * Child Patient Apron fits children under six years of age.

CODE	APRON DESCRIPTION
APA-3	Adult Patient Apron .30mm LE
APA-3-IT	Adult Patient Apron with Thyroid Shield .30mm LE
CPA-3	Child Patient Apron .30mm LE
CPA-3-IT	Child Patient Apron with Thyroid Shield .30mm LE
APA-5	Adult Patient Apron .50mm LE
APA-5-IT	Adult Patient Apron with Thyroid Shield .50mm LE
CPA-5	Child Patient Apron .50mm LE
CPA-5-IT	Child Patient Apron with Thyroid Shield .50mm LE



Pano Cape Apron

- * Drapes over shoulders.
- * Fastens in back with Velcro closure.
- * Protects shoulders and upper torso, front and back.

CODE	APRON DESCRIPTION
PC-3	Pano cape apron .30mm LE
PC-5	Pano cape apron .50mm LE



Pano Dual Apron (Adult/Child)

- * Two-panel.
- * Protects front and back from shoulders to gonadal areas.
- * Right shoulder Velcro closure.
- * Optional Velcro side closures

CODE	APRON DESCRIPTION
APD-3	Adult Pano Dual Apron .30mm LE
CPD-3	Child Pano Dual Apron .30mm LE
APD-5	Adult Pano Dual Apron .50mm LE
CPD-5	Child Pano Dual Apron .50mm LE