


# PATIENT LIFT

7600



**TECHLEM**  
STRETCHERS

SIMPLE AND EFFICIENT IN DESIGN, OUR PATIENT LIFT AND TRANSFER STRETCHER IS ADAPTABLE TO YOUR CLINICAL ENVIRONMENT. HANDCRAFTED FROM HIGH GRADE MATERIALS, THE 7600 IS HIGHLY DURABLE, MAINTAINING ITS FORM AND FUNCTION OVER DECADES OF USE.



HIGH DENSITY PATIENT TRANSFER SHEET, HIGH GRADE STEEL SHEET BUCKLES, AND TELESCOPING I.V. POLE INCLUDED

## MOVEMENT

HEIGHT ADJUSTMENT: 30" TO 42"  
5" ALL-SURFACE, LOCKING CASTERS  
5TH WHEEL FOR STEER ASSIST  
CRANK HANDLE FOR EASE OF USE

## DIMENSIONS

OUTSIDE DIMENSIONS: 77" X 31"  
INSIDE DIMENSIONS: 74" X 24<sup>1/2</sup>"

\* Weight Capacity: 300lbs