Apron Order Form

1.877.338.6857




Quote#:	User Name:		Date:
Rep.:	Facility Name:		Dept.:
APRON SPECS. Apron Code:		Qty:	Gender: Male (Apron Sizing, not User Gender) Female
Lead Type: [True Lead (150KVI	,	Cirrus ECO Lite (100KVP)
Color: Fab		FU-DF FU-PP Trim Inside	
Size: Coat Apron XS M2 S L M XL M1 2XL	Frontal Apron XS L S XL M 2XL	2-pc. Vest XS XL S 2XL M 3XL L	2-pc. Skirt XS M2 XL S L 2XL M1 L1 3XL M1 L2
APRON OPTIONS			
EMBROIDERY Script BLOCK Athletic	Embroidery: Location: Thread Colour: Name Tag:	1st Line Chest Skirt Thyroid Shield Other	2nd Line er Chest Skirt Thyroid Shield Other
	Thread Colour:		
THYROID SHIELD Standard BIB Visor Mammo Sewn-on(Integrated)	/ ATTACHMENT Loose Tethered Buckle Velcro Metal Ring Plastic Tie	Chest Dosimeter Skirt/Other	RM SHIELDS User Left User Right Sewn-on Snap-on Snap-on 10"/11"/12" Custom Size
Shipping Info/ Apron Details/ Modifications: Signed by:		Date	

The details of your apron order are recorded on the apron order form and these details are confirmed via your signature. Once confirmed, all subsequent changes are subject to charges. All aprons are customized in terms of style, fabric, embroidery and options and as such are non-returnable.

Made-To-Measure Apron Order Form





Instruction	ns:
	nstruction

-	

Vest Length	"+ 3" = total "		
Skirt Length	u		
Apron Length	ıı .		

Shoulder Width			"
Chest			"
Waist			"
Hips			"
Overall Height	1	"	

- Shoulders: Standing behind the person, measure the distance between the right and left clavicle bones. Do not measure end-of-shoulder to end-of-shoulder.
- Chest: Measure just under the arms and around the chest at the fullest point.
- Waist: Measure around the waistline where the apron waist will rest.
- Hip Size: Measure around the fullest point of the seat while standing.
- Overall Height : In addition to the above measurements,
 please provide the user's height to ensure a proper fit.

- Length (Coat & Frontal): Measure length from midshoulder to the knee (mid-knee or above-knee or belowknee) as per the user's preference or facility specifications.
- Length (Vest): Measure length from mid-shoulder to waist. Add an additional 3" for overlap.
- Length (Skirt) : Measure length from waist (where the apron waist will rest) to the knee (mid-knee or aboveknee or below-knee) as per the user's preference or facility specifications.