



ULTRARAY
MEDICAL


DENTAL: Radiation Protection and Accessories

Learn more at ultraraymedical.com

- ⊗ RADIOLOGY
- ⊗ CT and MRI
- ⊗ ELECTROMEDICAL
- ⊗ MAMMOGRAPHY
- ⊗ NUCLEAR MEDICINE
- ⊗ ULTRASOUND
- ⊗ ONCOLOGY/RADIATION THERAPY
- ⊗ STERILE PROCESSING
- ⊗ OR SUPPLIES
- ⊗ PATIENT HANDLING
- ⊗ FURNITURE

Embroidery Options

- Custom embroidery and stock designs available
- A logo or design specific to your facility will further enhance the product and clearly identify it.
- For more options, please refer to our main Medical Catalogue.



Thyroid Shields

- .30mm (for dental application only) or .50mm L.E.
- Available in all fabric colours.
- Can be ordered Loose or Integrated (sewn on).
- For integrated shields, please see apron styles.

CODE (0.30MM)	CODE (0.50MM)	DESCRIPTION
TSRL-3	TSRL-5	Standard Thyroid Shield Loose
TSML-3	TSML-5	Mammography STyle Thyroid Shield Loose



Technician Apron

- Full length: Protects from chest to knees
- 1" straps cross the back and tie securely at the front
- Available in XS-XXL, with or without integrated (sewn-on) thyroid collar.
- Easy on and off convenience.
- Optional no-charge chest pocket.
- Available in nylon or vinyl fabric options.



Tie-strap back



SIZE	SHOULDER WIDTH	MAX. CHEST/HIP	APRON LENGTH
X-Small	14"	34"	34"
Small	14.5"	38"	35"
Medium	16"	42"	37"
Large	18.5"	44"	39"
X-Large	20.5"	46"	41"
XX-Large	21.5"	50"	43"

CODE	DESCRIPTION
TA-3	Technician apron tie strap .30mm LE
TA-5	Technician apron tie strap .50mm LE

Patient Apron (Adult/Child)

- Available with or without thyroid collar.
- Hanging loops allow for proper storage.
- Child Patient Apron fits children under six years of age.

CODE	APRON DESCRIPTION
APA-3	Adult Patient Apron .30mm LE
APA-3-IT	Adult Patient Apron with Thyroid Shield .30mm LE
CPA-3	Child Patient Apron .30mm LE
CPA-3-IT	Child Patient Apron with Thyroid Shield .30mm LE
APA-5	Adult Patient Apron .50mm LE
APA-5-IT	Adult Patient Apron with Thyroid Shield .50mm LE
CPA-5	Child Patient Apron .50mm LE
CPA-5-IT	Child Patient Apron with Thyroid Shield .50mm LE



Pano Cape Apron

- Drapes over shoulders.
- Fastens in back with Velcro closure.
- Protects shoulders and upper torso, front and back.

CODE	APRON DESCRIPTION
PC-3	Pano cape apron .30mm LE
CPC-3	Child pano cape apron .30mm LE



Pano Dual Apron

- Two-panel
- Protects front and back from shoulders to gonadal areas.
- Right shoulder Velcro closure.
- Optional Velcro side closures.

CODE	APRON DESCRIPTION
APD-3	Adult pano dual apron .30mm LE
CPD-3	Child pano dual apron .30mm LE
APD-5	Adult pano dual apron .50mm LE
CPD-5	Child pano dual apron .50mm LE





1.877.338.6857 or info@ultraray.com
ultraraymedical.com
