



# Product Catalog

Patient Positioning

Diagnostic Imaging

Interventional Imaging

Operating Room Accessories

Distributed by



(877) 338-6857  
info@ultraray.com  
ultraraymedical.com

# Introduction

Domico Med-Device is a leading supplier of patient positioning products, accessories, and system components.

An effectively positioned patient reduces the risk of injury, minimizes the need for image retake and optimizes the efficiency of the modality. In this catalog you will be introduced to our patient positioning products with a wide variety of capabilities. With our team of highly-skilled engineers, designers, technicians and support professionals, Domico Med-Device fulfills customer needs while bringing passion to our pursuit of improving lives, together.

Questions? Call us at 810.750.5300 or email [orders@domicomed.com](mailto:orders@domicomed.com).

In addition to hundreds of procedure-specific and general use positioning devices, our standard foam positioning devices are available with the following alternatives:

## Standard Foam

- Open-cell foam
- Radiolucent
- Fire-retardant
- Economic general positioning

## Ultra Foam (-UF)

- Firm closed-cell foam
- Resists fluid absorption/contamination
- Fire-retardant
- Certain parts may require gluing to obtain complex dimensions (Glue joints may attenuate)

## Standard Vinyl Cover (-8)

- Frost-colored vinyl cover for standard foam positioning devices
- Heat-sealed to prevent fluid absorption/contamination (Seams may show in image if in field of view. Not recommended for low KVP techniques.)

## Heavy Duty Vinyl Cover (-20)

- Extremely durable vinyl cover for standard foam positioning devices
- Heat-sealed to prevent fluid absorption/contamination (Seams may show in image if in field of view. Not recommended for low KVP techniques.)

## ScanCoat® (-SCB – black)

- Strong, state-of-the-art sprayed vinyl coating:
- Withstands normal clinical use without tearing or peeling
- Very low attenuation
- Skid resistant
- Resists contamination (Edges may show in low KVP techniques)

Our product modalities include **General Positioning** , **CT Scan** , **MRI** , **Nuclear Medicine** , and **Ultrasound** . The icons shown here are used to represent each modality throughout this catalog.



## Medical Imaging

### Featured Products

Spine Tabletop .....	4
Aquilion CT Accessories .....	5

### Positioning Aids

Torso/Body .....	6-7
Finger .....	7
Leg .....	8
Head and Neck .....	9-10
Positioning Kits .....	11-14
Pediatric Positioning .....	14-15
Rectangles & Wedges .....	16-18

Table Pads .....	19
------------------	----

Specialty Procedure Items .....	20
---------------------------------	----

Armboards .....	20-21
-----------------	-------

MRI Accessories .....	21
-----------------------	----

MEDVAC Immobilization System .....	22
------------------------------------	----

CT Accessories .....	22
----------------------	----

### Canon Aquilion CT Accessories

Short System .....	23
Long System .....	24
One System .....	25
Accessories for All Systems .....	26
Head Holder Kit .....	26-27
Table Pads & Covers .....	28-29
Detachable Side Rails & Accessories .....	30-32

Carbon Fiber Head Holder .....	32
--------------------------------	----

Veterinary Positioning Aids .....	33-34
-----------------------------------	-------

## O.R. Accessories

Urology Accessories .....	35
---------------------------	----


Not all products are in stock and lead times may vary.

# Spine Tabletop

**Part Number 32588**  
**Spine Tabletop with accessories**

Compatible with Siemens Artis Zeego,  
Artis Q, Artis Zee

**Procedures** • Lumbar Fusion • SI Fusion • Thoracic Fusion • Lumbar Decompression • Anterior Cervical Discectomy and Fusion • Posterior Cervical Fusion • Extreme Lateral Interbody Fusion • Oblique Lateral Interbody Fusion • Brain Biopsy • Ortho Trauma



- 1 Radiolucent Carbon Fiber Construction**   **2 Spar configuration allows the patient's abdomen to hang free**
- 3 Three pad sizes (S, M, L) allow maximum options for patient and surgeon comfort**   **4 Dovetail adapter allows for Doro Head Clamp integration**
- 5 Face Plate mirror allows the surgeon to monitor patient's alertness/awareness**
- 6 Disposable face pad provides single patient use comfort**   **7 Straps secure patient to tabletop**


## Increase Interventional X-Ray Utilization

Adopting and integrating the Spine Tabletop into your hybrid OR suite will increase capacity and efficiency by providing spinal and neurosurgeons a new venue to perform minimally invasive procedures. Install in 4 easy steps.

- Designed for minimally invasive spine procedures performed with the patient prone.
- No new software programming required. Tabletop functions using the existing Neuro Tabletop program.
- Removable design allows simple attachment to the pedestal for fast room conversion.



# Aquilion CT Accessories










See complete product listings on pages 23–27.











**NOTE:** Gray security straps, with male fitting which couples to a female security rail, are not Domico Med-Device products. If customer is using gray security straps and orders new straps from Domico Med-Device (9061, 9062) they must order the corresponding detachable rail w/extender for their table. 9078, 9079 or 9080.

# Torso/Body Positioning



















Product Number	Product Description	Application	Alternatives	Modality	
100	<b>Torso Block 27°</b> 7"H x 16"W x 12"L	Dorsal and lumbar spine support	-UF, -8, -20, -SCB	    	 100
107	<b>Angular Bolster</b> 7.25"H x 12.5"W x 21"L	Torso and aortogram procedures	-UF, -8, -20, -SCB	    	 107
109	<b>Angular Bolster</b> 7.25"H x 12.5"W x 15"L	Mammograph procedures	-UF, -8, -20, -SCB	    	 109
115	<b>Bolster</b> 5"H x 7"W x 15"L	Knee and back support	-UF, -8, -20, -SCB	    	 115
124	<b>Bolster</b> 6.6"H x 9"W x 18"L	Stabilization for myelogram needle insertion	-UF, -8, -20, -SCB	    	 124
145	<b>Spinal Wedge</b> 7"H x 7"W x 21.25"L	Oblique spine support	-UF, -8, -20, -SCB	    	 145
146	<b>Cylinder</b> 6"D x 12"L	Abdominal compressions during kidney biopsies	-UF, -8, -20, -SCB	    	 146
153	<b>Cylinder</b> 4.25"D x 12"L	Abdominal compression during kidney biopsies	-UF, -8, -20, -SCB	    	 153
150	<b>Myelography Block</b>	Abdomen elevation for proper contrast flow	-SCB		 150
160	<b>Lumbar Support</b> 3"H x 10"W x 10"L	Spine positioning	-UF, -8, -20, -SCB	    	 160-161-162
161	<b>Lumbar Support</b> 4"H x 10"W x 10"L	Spine positioning	-UF, -8, -20, -SCB	    	
162	<b>Lumbar Support</b> 5"H x 10"W x 10"L	Spine positioning	-UF, -8, -20, -SCB	    	

Product Number	Product Description	Application	Alternatives	Modality	
195	<b>Cardiac Wedge</b> 14"H x 28"W x 16"L	Cardiac procedures	-UF, -8, -20, -SCB	  	
438	<b>Rounded Shoulder Support</b> 9.5"H x 17"W x 20.75"L	Provides comfort when patient must extend arms out of the field of view while lying supine	-SCB		 438
439	<b>Winged Shoulder Support</b> 10.3"H x 19"W x 14.6"L	Provides comfort when patient must extend arms out of the field of view while lying supine	-SCB		 439

## Hand Positioning

Product Number	Product Description	Application	Alternatives	Modality	
121	<b>Oblique Finger Block</b>	Right and left oblique hand imaging	-UF, -SCB		 121
428	<b>Lateral Finger Block Right Hand</b>	Positions fingers and wrist in a lateral plane. Excellent for individual finger studies.	-UF, -SCB		 428
429	<b>Lateral Finger Block Left Hand</b>	Positions fingers and wrist in a lateral plane. Excellent for individual finger studies.	-UF, -SCB		 429
430	<b>Dual Oblique Finger Block</b>	Right and left oblique hand imaging	-UF, -SCB		 430
431	<b>Dual Lateral Finger Block</b>	Positions finger and wrist in a lateral plane. Excellent for individual finger studies on both hands simultaneously.	-UF, -SCB		 431

# Leg Positioning
















Product Number	Product Description	Application	Alternatives	Modality	
203	<b>Cross-Table Lateral Leg Block</b>	Used for comparative studies of legs	-UF, -SCB		 203
209	<b>Knee Arthrography Block</b> 8"H x 12"W x 10"L	Facilitates immobilization and effective knee flexion during arthrography.	-UF, -8, -20, -SCB		 209
401	<b>Leg Immobilizer</b>	Positions and immobilizes lower extremities during venogram or arteriogram exams.	-UF, -SCB	    	 401
427	<b>Leg Brace</b> 7" x 10" x 20"	Contoured cavity keeps leg in place. Used with a cast for comfort.	-UF, -SCB	    	 427
466	<b>Sloped Leg Immobilizer</b>	Used for arthogram exams	-SCB		 466










# Head and Neck Positioning

Product Number	Product Description	Application	Alternatives	Modality	
101	<b>Basal Block</b> 7"H x 16"W x 12"L	Skull work as well a lateral knees, ankles and legs	-UF, -SCB	   	 101
102	<b>AP-PA Head Holder</b>	Holds patient's head in several positions allowing multiple shots	-UF, -SCB	   	 102
104	<b>Pressure Block</b> 5.5"H x 7.3"W x 10"L	Special procedure skull work	-UF, -8, -20, -SCB	    	 104
112	<b>Circular Disc</b> 1.5"H x 7" Diameter	Lifts the skull to properly position the spine	-UF, -8, -20, -SCB	    	 112
149	<b>Circular Disc</b> 2.5"H x 7" Diameter	Lifts the skull to properly position the spine	-UF, -8, -20, -SCB	    	 149
163	<b>Concave Disc</b> 1.5"H x 7" Diameter	General skull work	-SCB	    	 163-164-165
164	<b>Concave Disc</b> 2.5"H x 7" Diameter	General skull work	-SCB	    	
165	<b>Concave Disc</b> 3"H x 7" Diameter	General skull work	-SCB	    	
119	<b>Head and Neck Rest</b>	Special procedure skull work. Also therapy and emergency situations.	-UF, -SCB	    	 119
132	<b>Head and Neck Support</b>	General skull work	-UF, -SCB	    	 132
133	<b>Cervical Head Rest</b> 4.25"H x 4.75"W x 8"L	General skull work	-UF, -8, -20, -SCB	    	 133
205	<b>Oblique Cervical Skull Block</b>	Postions skull comfortably at 45° oblique spine	-UF, -SCB		 205
207	<b>4-in-1 Lateral Skull Cube</b>	Flip cube to appropriate height 5", 6" 7" & 8"	-SCB		 207

# MEDICAL IMAGING POSITIONING AIDS

Product Number	Product Description	Application	Alternatives	Modality	
405	<b>Waters Positioner</b>	Laminographs of the orbits and facial bones	-UF, -SCB	 	 405
407	<b>Dual Lateral Skull Block</b>	Provides both right and left lateral views. Used for trauma, routine skull work and laminographs	-SCB		 407
408	<b>AP-PA Head Immobilizer</b>	Immobilize traumatized patients. Excellent for ER use, can be taped to a table.	-UF, -SCB	   	 408
409	<b>AP Head Immobilizer</b>	Includes a hook & loop strap. Used for positioning or immobilization.	-SCB	   	 409

## Head and Neck Misc.

Product Number	Product Description	Application	Alternatives	Modality	
300	<b>Head Paddles</b> Packaged in Pairs	Thick foam pads rigid plastic bases. Designed with oversized handles for easy, secure grip.	-8, -SCB		 300
303	<b>Head-Ends</b> Packaged in Pairs	Unique skull immobilization concept. Fifteen pound weights are fastened to vinyl covered foam pads. Non-skid rubber pads on the bottom prevent them from moving.	-UF		 303
9025	<b>Coronal Head Holder</b>	Allows patients to rest comfortably in prone position. Made of clear, 1/4" acrylic with two coated foam pads.		 	 9025

# Positioning Kits

Product Number	Product Description	Includes	Alternatives	Modality	
1	<b>Deluxe Kit</b> 13 pieces 9 part numbers	(2) 110 - wedges (1) 111 - wedge (1) 112 - circular disc (2) 113 - wedges (1) 114 - rectangle (1) 151 - decubitus block (2) 153 - cylinder (2) 154 - thin mattresses (1) 171 - table pad	-UF, -8, -20, - SCB		 Kit 1
2	<b>Large Kit</b> 8 pieces 7 part numbers	(1) 100 - torso block (1) 101 - basalar blk (1) 102 - adult AP-PA head holder (1) 104 - pressure block (2) 105 - wedges (1) 106 - wedge (1) 169 - rectangle	-UF, - SCB		 Kit 2
3	<b>Small Kit</b> 7 pieces 5 part numbers	(2) 110 - wedges (1) 111 - wedge (1) 112 - circular disc (2) 113 - wedges (1) 114 - rectangle	-UF, -8, -20, - SCB		 Kit 3
4	<b>Kit A</b> 7 pieces 5 part numbers	(1) 112 - circular disc (1) 114 - rectangle (2) 116 - wedges (2) 117 - wedges (1) 118 - wedge	-UF, -8, -20, - SCB		 Kit 4
5	<b>Kit B</b> 8 pieces 6 part numbers	(1) 112 - circular disc (1) 114 - rectangle (2) 116 - wedges (2) 117 - wedges (1) 118 - wedge (1) 121 - oblique finger block	-UF, - SCB		 Kit 5
6	<b>Kit C</b> 5 pieces 5 part numbers	(1) 100 - torso block (1) 102 - adult AP-PA head holder (1) 104 - pressure block (1) 114 - rectangle (1) 121 - oblique finger block	-UF, - SCB		 Kit 6

# MEDICAL IMAGING POSITIONING AIDS

Product Number	Product Description	Includes	Alternatives	Modality	
7	<b>Kit D</b> 4 pieces 4 part numbers	(1) 128 - wedge (1) 129 - wedge (1) 130 - wedge (1) 140 - square	-UF, -8, -20, - SCB		 Kit 7
8	<b>Kit E</b> 10 pieces 10 part numbers	(1) 101 - basalar block (1) 102 - adult AP-PA head holder (1) 104 - pressure block (1) 128 - wedge (1) 129 - wedge (1) 130 - wedge (1) 131 - square (1) 136 - wedge (1) 138 - wedge (1) 140 - wedge	-UF, - SCB		 Kit 8
9	<b>Rectangle Kit</b> 12 pieces	(3) .5" x 4" x 6" (3) 1" x 4" x 6" (3) 1.5" x 4" x 6" (3) 2" x 4" x 6"	-UF, -8, -20, - SCB		 Kit 9
10	<b>Kit F</b> 6 pieces 6 part numbers	(1) 101 - basalar block (1) 102 - adult AP-PA head holder (1) 104 - pressure block (1) 131 - square (1) 136 - wedge (1) 138 - wedge	-UF, - SCB		 Kit 10
11	<b>Special Procedure Skull Kit</b> 4 pieces 4 part numbers	(1) 405 - waters positioner (1) 406 - ear laminograph block (1) 407 - dual lateral skull block (1) 408 - AP-PA head immobilizer	-UF, - SCB		 Kit 11
14	<b>Ortho Kit</b> 5 pieces 5 part numbers	(1) 102 - adult AP-PA head holder (1) 121 - oblique finger block (1) 133 - cervical head rest (1) 145 - wedge (1) 169 - rectangle	-UF, - SCB		 Kit 14

Product Number	Product Description	Includes	Alternatives	Modality	
15	<b>Kit I</b> 7 pieces 7 part numbers	(1) 127 - child AP-PA head holder (1) 132 - head/neck support (1) 133 - cervical head rest (1) 136 - wedge (1) 138 - wedge (1) 139 - wedge (1) 169 - rectangle	-UF, - SCB		 Kit 15
17	<b>Routine X-Ray Kit</b> 10 pieces 7 part numbers	(1) 112 - circular disc (1) 114 - rectangle (2) 116 - wedges (2) 117 - wedges (1) 118 - wedge (1) 121 - oblique finger block (2) 122 - wedges	-UF, - SCB		 Kit 17
19	<b>Standard R &amp; F Kit</b> 8 pieces 7 part numbers	(1) 111 - wedge (1) 112 - circular disc (2) 113 - wedges (1) 114 - rectangle (1) 121 - oblique finger block (1) 133 - cervical head rest (1) 145 - wedge	-UF, - SCB		 Kit 19
99	<b>CT Kit</b> 10 pieces 8 part numbers	(2) 105 - wedges (1) 106 - wedge (2) 113 - wedges (1) 118 - wedge (1) 120 - wedge (1) 123 - wedge (1) 130 - wedge (1) 136 - wedge (1) 138 - wedge (1) 140 - wedge	-UF, -8, -20, - SCB		 Kit 99

# Positioning Kit Accessories

Product Number	Product Description	Application
422	<b>Block Hanger Kit</b>	Provides fast on/off the wall storage of Position Perfect positioning aids. The kit includes Velstick and 8 Velcoins for adhesion.
423	<b>Velcoins</b> set of 24 13 pieces	24 adhesive Velcro coins on a convenient peel off sheet.



422



423

# Pediatric Positioning



































































Product Number	Product Description	Application	Alternatives	Modality		
VMR33PUBB	<b>Pediatric Immobilizer</b> Newborn to 1 year old	Safely and securely fastens around the patient restricting movement and reducing repeat scans.				VMR33PUBB, VMR38X01
VMR38X01	<b>Child Immobilizer</b> 1 year to 4 years old	Safely and securely fastens around the patient restricting movement and reducing repeat scans.				
425	<b>The Pediatric Hugger Complete Unit</b> Up to 3 months	Provides womb-like warmth and security during immobilization.				425-426
426	<b>The Child Hugger Complete Unit</b> 4 months to 3 years	Provides womb-like warmth and security during immobilization.				
200	<b>CT Hugger Complete Unit</b>	Radiused base prevents artifacts. Aluminum equivalency below 1.5mm.				200
127	<b>Child AP-PA Head Holder</b>	Comfortably holds head in multiple positions allowing multiple shots.	-UF, -8, -20, -SCB			127

Product Number	Product Description	Application	Alternatives	Modality	
148	<b>Child Myelography Block</b>	Elevates abdomen for proper contrast media flow.	-SCB		 148
301	<b>Pediatric Dual Oblique Finger Block</b>	Separation of fingers for overhead lateral views of both hands or individual finger studies.	-UF		 301
302	<b>Pediatric Dual Lateral Finger Block</b>	45° oblique of both hands for comparative studies or individual finger studies.	-UF		 302
304	<b>AP-PA Head Immobilizer - Pediatric</b>	Both pediatric and child AP-PA head immobilizers provide AP-PA views where extreme weight is not necessary. Excellent in emergency areas. May be Velcro fastened to table.	-UF, -SCB	  	 304
305	<b>AP-PA Head Immobilizer - Child</b>	Both pediatric and child AP-PA head immobilizers provide AP-PA views where extreme weight is not necessary. Excellent in emergency areas. May be Velcro fastened to table.	-UF, -SCB	  	 305
308	<b>Pediatric AP-PA Hand/Foot Immobilizer</b>	Immobilizes right and left hands or feet with pressure conforming foam for comparative studies. Soft, clear vinyl compression band provides correct positioning during examination. For use with 8" x 10" film cassette. Vinyl components are heat-sealed to prevent contamination and resist stains. Rubber pads on bottom of 6" x 10" Plexiglas base, provide stability.			 308

# Positioning Rectangles & Wedges

Product Number	Product Description	Application	Alternatives	Modality	
103	<b>Rectangle</b> 3"H x 8"W x 12"L	General Positioning	-UF, -8, -20, -SCB	    	 103
114	<b>Rectangle</b> 2"H x 7.5"W x 9.5"L	General Positioning	-UF, -8, -20, -SCB	    	 114
126	<b>Rectangle - Set of 4</b> (1) .5" x 4" x 6" (1) 1" x 4" x 6" (1) 1.5" x 4" x 6" (1) 2" x 4" x 6"	General Positioning	-UF, -8, -20, -SCB	    	 126
131	<b>Rectangle</b> 4"H x 10"W x 10"L	General Positioning	-UF, -8, -20, -SCB	    	 131
142	<b>Rectangle</b> 3"H x 4.75"W x 9.5"L	General Positioning	-UF, -8, -20, -SCB	    	 151
143	<b>Rectangle</b> 4.75"H x 4.75"W x 9.5"L	General Positioning	-UF, -8, -20, -SCB	    	 151
151	<b>Rectangle</b> 4"H x 18"W x 24"L	General Positioning	-UF, -8, -20, -SCB	    	 154
152	<b>Decubitus Block</b> 3"H x 14"W x 18"L	General Positioning	-UF, -8, -20, -SCB	    	 154
154	<b>Thin Mats - Set of 2</b> .25"H x 18"W x 24"L	General Positioning	-UF, -8	    	 169
169	<b>Decubitus Block</b> 3"H x 7.5"W x 10"L	General Positioning	-UF, -8, -20, -SCB	    	 169
105	<b>Wedge - 18.75°</b> 3.75"H x 9"W x 8"L	General Positioning	-UF, -8, -20, -SCB	    	 105
106	<b>Wedge - 45°</b> 6.6"H x 6.6"W x 12"L	General Positioning	-UF, -8, -20, -SCB	    	 106
108	<b>Wedge - 45°</b> 5"H x 5"W x 15"L	General Positioning	-UF, -8, -20, -SCB	    	 106



Product Number	Product Description	Application	Alternatives	Modality	
110	<b>Wedge - 45°</b> Set of 2 3"H x 3"W x 7.25"L	General Positioning	-UF, -8, -20, -SCB	    	 110
111	<b>Wedge - 45°</b> 6.75"H x 6.75"W x 7.25"L	General Positioning	-UF, -8, -20, -SCB	    	 111
113	<b>Wedge - 12°</b> Set of 2 2.25"H x 9.5"W x 7.25"L	General Positioning	-UF, -8, -20, -SCB	    	 113
116	<b>Wedge - 34°</b> Set of 2 3"H x 4"W x 7"L	General Positioning	-UF, -8, -20, -SCB	    	 116
117	<b>Wedge - 16°</b> Set of 2 2.5"H x 7"W x 10"L	General Positioning	-UF, -8, -20, -SCB	    	 117
118	<b>Wedge - 45°</b> 7"H x 7"W x 7"L	General Positioning	-UF, -8, -20, -SCB	    	 118
120	<b>Wedge - 18°</b> 2.5"H x 8"W x 7"L	General Positioning	-UF, -8, -20, -SCB	    	 120
122	<b>Wedge - 22°</b> Set of 2 3.5"H x 7"W x 8"L	General Positioning	-UF, -8, -20, -SCB	    	 122
123	<b>Wedge - 33°</b> 6.6"H x 9"W x 12"L	General Positioning	-UF, -8, -20, -SCB	    	 123
128	<b>Wedge - 13.5°</b> 3"H x 10"W x 10"L	General Positioning	-UF, -8, -20, -SCB	    	 128
129	<b>Wedge - 21°</b> 4.6"H x 10"W x 10"L	General Positioning	-UF, -8, -20, -SCB	    	 129
130	<b>Wedge - 45°</b> 10"H x 10"W x 10"L	General Positioning	-UF, -8, -20, -SCB	    	 130
136	<b>Wedge - 12.5°</b> 2.25"H x 8"W x 10"L	General Positioning	-UF, -8, -20, -SCB	    	 136
138	<b>Wedge - 21°</b> 3.75"H x 8"W x 10"L	General Positioning	-UF, -8, -20, -SCB	    	 138

# MEDICAL IMAGING POSITIONING AIDS

Product Number	Product Description	Application	Alternatives	Modality	
139	<b>Wedge - 45°</b> 7"H x 7"W x 10"L	General Positioning	-UF, -8, -20, -SCB	    	
140	<b>Wedge - 27°</b> 6"H x 10"W x 10"L	General Positioning	-UF, -8, -20, -SCB	    	
141	<b>Spinal Wedge - 45°</b> 9"H x 9"W x 30"L	General Positioning	-UF, -8, -20, -SCB	    	141
144	<b>Wedge - 45°</b> 7"H x 7"W x 14"L	General Positioning	-UF, -8, -20, -SCB	    	
145	<b>Spinal Wedge - 45°</b> 7"H x 7"W x 21.25"L	General Positioning	-UF, -8, -20, -SCB	    	144
158	<b>Wedge - 15°</b> <b>Set of 2</b> 3"H x 11"W x 7"L	General Positioning	-UF, -8, -20, -SCB	    	
159	<b>Wedge - 20°</b> <b>Set of 2</b> 4"H x 11"W x 7"L	General Positioning	-UF, -8, -20, -SCB	    	145
191	<b>Wedge - 45°</b> 7"H x 7"W x 24"L	General Positioning	-UF, -8, -20, -SCB	    	
196	<b>Wedge - 35/55°</b> 9.75" x 6.75" x 24"L	General Positioning	-UF, -8, -20, -SCB	    	
201	<b>Wedge - 35/55°</b> 11"H x 7"W x 30"L	General Positioning	-UF, -8, -20, -SCB	    	191









## Standard Table Pads & Covers

Product Number	Product Description	Product Description	Size
Foam Pad Only	Foam Pad with -8 Cover	Foam Pad with -20 Heavy Duty Cover	
170	170-180	170-20	1.5" x 21" x 42"
171	171-181	171-20	1" x 24" x 72"
172	172-182	172-20	1.25" x 26" x 78"
173	173-183	173-20	1.5" x 27" x 78"
174	174-184	174-20	1" x 21" x 42"
175	175-185	175-20	1" x 24" x 60"
176	176-186	176-20	1" x 27" x 78"



171-176







## Premium Heavy Duty Table Pads

Product Number	Product Description	Application	Modality
520	Table Pad with heavy duty 20 mil vinyl cover 1" x 24" x 72"	Grommets are located at the top for hanging. Cross-hair centering helps to easily position the patient.	   
530	Table Pad with heavy duty 20 mil vinyl cover 1" x 19" x 72"	Grommets are located at the top for hanging. Cross-hair centering helps to easily position the patient.	   





520-530

# Specialty Procedure Items

Product Number	Product Description	Application	Alternatives	Modality	
135	<b>EZ Straps</b> 1.5" x 30' roll White hook & loop	Cut to desired length.			 135
168	<b>Security Strap</b> 5"W x 72"L	Secures patient to table or wheelchair.			 168
199	<b>Weight-Bearing Cassette Holder</b>	Facilitates longitudinal arch and lateral projection studies.			 199
204	<b>Decubitus Block w/ Film Cassette Holder</b>	Support and spine correction for laterals and decubitus views.			 204
208	<b>Auto-Traction Straps</b> 50"L with 14" Adjustable Wrist Straps	Images cervical vertebrae, in a lateral view, supine position. Patient applies equal traction on each arm by extending legs, thus tightening the straps.			 208
215	<b>Endovaginal Wedge</b> 6"H x 20"W x 20"L	Improves transducer access into vaginal canal. Heavy duty, fluid resistant cover is easily cleaned.			 215

## Armboards

Product Number	Product Description	Application	
250	<b>Carbon Fiber Armboard Rail Mount</b>	Properly positions patient. Rotate to desired position, helping to ease radial and brachial access. Optional pad, part 258, available	 250
258	<b>Armboard Pad Fits 250 &amp; 254</b>	Fits over armboards 250 & 254. Made with memory foam for extra comfort.	 258

Product Number	Product Description	Application
254	<b>Carbon Fiber Armboard Shoulder Mount</b>	Fits between table and table pad. Properly positions patient for radial and brachial access. Optional pad, part 254, available.
261	<b>Rest Eazy Armboard GE</b>	Fits GE tables. Positions patient's arms while conducting prone procedures on angio table.
262	<b>Rest Eazy Armboard Siemens &amp; Philips</b>	Fits Siemens & Philips tables. Positions patient's arms while conducting prone procedures on table.
263	<b>Rest Eazy Armboard Canon</b>	Fits most Canon tables. Positions patient's arms while conducting prone procedures on angio table.
265	<b>Rest Eazy Armboard Canon 880B Table</b>	Fits Canon 880B tables. Positions patient's arms while conducting prone procedures on angio table.



## MRI Accessories

Product Number	Product Description	Application	Alternatives
463	<b>MRI Comfort Pillow</b>	Comfortably supports the patient's head during procedures where the arms must be overhead.	
192	<b>Leg Rest with Knee Cutouts</b> Small 5" x 18" x 15"	Provides leg and knee support while maintaining hip flexion. Helps relieve lower back stress and pain.	-SCB
193	<b>Leg Rest with Knee Cutouts</b> Medium 7" x 13" x 17.5"	Provides leg and knee support while maintaining hip flexion. Helps relieve lower back stress and pain.	-SCB
194	<b>Leg Rest with Knee</b> Large 10" x 24" x 18"	Provides leg and knee support while maintaining hip flexion. Helps relieve lower back stress and pain.	-SCB



# MedVac Immobilization System

Product Number	Product Description	Application
VMR33PUBB	Pediatric Splint - Newborn to 1 year old	
VMR38X01	Child Splint - 1 year to 4 years old	
27625	Manual Pump	



## CT Accessories

Product Number	Product Description	Application
99	<p><b>CT Kit</b></p> <p>(2) 105 - wedge 3.75" H x 9" W x 8" L                      (1) 106 - wedge 6.6" H x 6.6" W x 12" L                      (2) 113 - wedge 2.25" H x 9.5" W x 7.25" L                      (1) 118 - wedge 7" H x 7" W x 7" L                      (1) 120 - wedge 2.5" H x 8" W x 7" L                      (1) 123 - wedge 6.6" H x 9" W x 12" L                      (1) 130 - wedge 10" H x 10" W x 10" L                      (1) 136 - wedge 2.25" H x 8" W x 10" L                      (1) 138 - wedge 3.75" H x 8" W x 10" L                      (1) 140 - wedge 6" H x 10" W x 10" L</p>	



**156 CT Pillow**  
3" H x 10" W x 12" L Tapered edges will not affect imaging.

**192 Leg Rest with Knee Cutouts**  
Small 5" x 18" x 15" Provides leg and knee support while maintaining hip flexion. Helps relieve lower back stress and pain.

**193 Leg Rest with Knee Cutouts**  
Medium 7" x 13" x 17.5" Provides leg and knee support while maintaining hip flexion. Helps relieve lower back stress and pain.

**194 Leg Rest with Knee Cutouts**  
Large 10" x 24" x 18" Provides leg and knee support while maintaining hip flexion. Helps relieve lower back stress and pain.

**200 CT Hugger Complete Unit**  
Provides patients newborn to 3 years old) with womb-like warmth and security during immobilization. CT and MRI safe. Radiused base helps prevent artifacts.



# Aquilion Accessories Short System

Product Number	Product Description	Application
9100	<b>Aquilion Standard Accessory Kit-Short with Head Holder</b>	For use with Aquilion short system. Includes: 9106F, 9109, 9078, 9090, 9103F
9106F	<b>Short Table Pad</b>	Table Pad for Aquilion short couch, 1.00 x 18.60 x 73.80 with scan range line
9109	<b>Short Table Pad Cover</b>	Frosted vinyl table cover "sleeve" with flap. Access holes on bottom side of cover allow pad to adhere to table top Velcro. Compatible with all previous table pad designs.*
9078	<b>Detachable Slide Rail with Extender Short Slide Rail</b>	Attaches to couch separately from pad. Companion to the slide straps 9061 and 9062 for easy strap positioning within scan range only. 57". Packed inside table pad carton.**
9090	<b>Aquilion Pediatric &amp; Adult Head Holder Kit</b>	Accommodate all patients with a single head holder. All patient ages are comfortably positioned for accurate imaging using any of three sizes of head holders and wedge pads placed behind and/or beside the head.
9103F	<b>Aquilion Standard Accessory Kit-Short without Head Holder</b>	Identical to 9100 - Aquilion Standard Accessory Kit-Short but NOT including 9090 - Aquilion. Pediatric & Adult Head Holder Kit



9106



9109



9078



9090

\* Previous versions of the table cover are no longer available. If your current table pad has flaps they must be removed to accommodate the new table pad by simply cutting them off. This will not damage the pad or imaging quality. \*\* Aquilion systems are often shipped, not only with Domico Med-Device straps but with additional gray strap sets with a male runner. These are not compatible with Domico Med-Device parts and cannot be re-ordered from Domico Med-Device.

# Aquilion Accessories Long System

Product Number	Product Description	Application
9101	<b>Aquilion Standard Accessory Kit-Long with Head Holder</b>	For use with Aquilion long system. Includes: 9107, 9110, 9079, 9090, 9104
9107	<b>Long Table Pad</b>	Table Pad for Aquilion long couch, 1.00 x 18.60 x 85.60 with scan range line
9110	<b>Long Table Pad Cover</b>	Frosted vinyl table cover with flap. Access holes on bottom side of cover allow pad to adhere to table top Velcro. Compatible with all previous table pad designs.*
9079	<b>Detachable Slide Rail with Extender Long</b>	Slide rail attaches to couch separately from pad. Companion to the slide straps 9061 and 9062 for easy strap positioning within scan range only. 69". Packed inside table pad carton.**
9090	<b>Aquilion Pediatric &amp; Adult Head Holder Kit</b>	Accommodate all patients with a single head holder. All patient ages are comfortably positioned for accurate imaging using any of three sizes of head holders and wedge pads placed behind and/or beside the head.
9104	<b>Aquilion Standard Accessory Kit-Short without Head Holder</b>	Identical to 9101 - Aquilion Standard Accessory Kit-Short but NOT including 9090 - Aquilion. Pediatric & Adult Head Holder Kit.



9107



9110



9079



9090

\* Previous versions of the table cover are no longer available. If your current table pad has flaps they must be removed to accommodate the new table pad by simply cutting them off. This will not damage the pad or imaging quality. \*\* Aquilion systems are often shipped, not only with Domico Med-Device straps but with additional gray strap sets with a male runner. These are not compatible with Domico Med-Device parts and cannot be re-ordered from Domico Med-Device.




# Aquilion Accessories One System

Product Number	Product Description	Application
9102	<b>Aquilion High Capacity Couch Accessory Kit with Head Holder</b>	For use with Aquilion one system. Includes: 9108, 9111, 9080, 9090, 9105F
9108	<b>One Table Pad</b>	Table Pad for Aquilio One 1.00 x 18.60 x 96.30 with scan range line
9111	<b>One Table Pad Cover</b>	Frosted vinyl table cover with flap. Access holes on bottom side of cover allow pad to adhere to table top Velcro. Compatible with all previous table pad designs.*
9080	<b>Detachable Slide Rail with Extender One</b>	Slide rail attaches to couch separately from pad. Companion to the slide straps 9061 and 9062 for easy strap positioning within scan range only. 77". Packed inside table pad carton.**
9090	<b>Aquilion Pediatric &amp; Adult Head Holder Kit</b>	Accommodate all patients with a single head holder. All patient ages are comfortably positioned for accurate imaging using any of three sizes of head holders and wedge pads placed behind and/or beside the head.
9105F	<b>Aquilion High Capacity Couch Accessory Kit without Head Holder</b>	Identical to the above 9102 Aquilion High Capacity Couch Accessory Kit but NOT including 9090 Aquilion Pediatric & Adult Head Holder Kit



9108



9111



9080



9090

\* Previous versions of the table cover are no longer available. If your current table pad has flaps they must be removed to accommodate the new table pad by simply cutting them off. This will not damage the pad or imaging quality. \*\* Aquilion systems are often shipped, not only with Domico Med-Device straps but with additional gray strap sets with a male runner. These are not compatible with Domico Med-Device parts and cannot be re-ordered from Domico Med-Device.

# Aquilion Accessories Use With All Aquilion Systems

Product Number	Product Description	Application
9090	<b>Aquilion Pediatric &amp; Adult Head Holder Kit</b>	Accommodate all patients with a single head holder. All patient ages are comfortably positioned for accurate imaging using any of three sizes of head holders and wedge pads placed behind and/or beside the head. Includes: 9091, 9092, 9093, 9094, 9095, 9096, 9097
9091	<b>Carbon Fiber Head Holder</b>	For use with Aquilion Short, Long, or One. Head holder, carbon fiber: head pads and wedges accommodate all patients.
9092	<b>Adult Head Pad</b>	For use with Aquilion Short, Long, or One. Head pad, adult
9093	<b>Child Head Pad</b>	For use with Aquilion Short, Long, or One. Head pad, child
9094	<b>Pediatric Head Pad</b>	For use with Aquilion Short, Long, or One. Head pad, Pediatric
9095	<b>Side Wedge Pads</b>	For use with Aquilion Short, Long, or One. Wedge, side: for use between side of patient's head/cheek and with head pads
9096	<b>Bottom Wedge Pads</b>	For use with Aquilion Short, Long, or One. Wedge, bottom: for use with head pad and bottom of head holder, elevates patient's forehead
9097	<b>Strap Set of 3</b>	For use with Aquilion Short, Long, or One. Head straps, set of three. For use with 9090 and 9073 kit



9091



9092, 9093, 9094







## Aquilion Pediatric & Adult Head Holder Kit





Product Number	Product Description	Application
9090	<b>Aquilion Pediatric &amp; Adult Head Holder Kit</b>	Accommodate all patients with a single head holder. All patient ages are comfortably positioned for accurate imaging using any of three sizes of head holders wedge pads placed behind and/or beside the head.

9006	<b>Chin Strap</b>	For use with Aquilion Short, Long, or One. Works with Toshiba head holder and 9091, comfortably secures patient's head.	
9007	<b>Forehead Strap with Adult Pad</b>	For use with Aquilion Short, Long, or One. Ensure excellent imaging, patient comfort, and secure positioning using the 9090 and Toshiba head holders.	9006
9008	<b>Adult Head Rest Kit</b>	For use with Aquilion Short, Long, or One. These soft, image-friendly, pads fit inside the Toshiba head holders increasing patient comfort and stabilization.	
9009	<b>Tilt Wedge</b>	For use with Aquilion Short, Long, or One. Unique, angled, wedge designed for procedure-specific positioning needs.	
9017	<b>Knee Wedge</b>	For use with Aquilion Short, Long, or One. Coated, foam, wedge comfortably elevates patient's knees while maintaining image integrity.	9017
9020F	<b>Pediatric Lift Pad</b>	For use with Aquilion Short, Long, or One. Provides additional padding for children whose shoulders may feel uncomfortable on the head holder attachment points.	
9025	<b>Coronal Head Extremity Holder</b>	For use with Aquilion Short, Long, or One. This comfortable, prone-imaging, head holder provides unobstructed imaging of the temples. Also serves as a foot positioning device by placing the heel in the chin pad and the ball of the foot on the forehead pad.	
9061	<b>Wide Security Strap</b>	For use with Aquilion Short, Long, or One. Black, nylon, straps glide along slide rails on each side of table; affix to one another for easy securing of patients. Two at 14" **	9025
9113	<b>Wide Bariatric Extension Strap</b>	For use with Aquilion Short, Long, or One. Affixes and lengthens 9061 Wide Security Strap for bariatric patients.**	
9062	<b>Medium Security Strap</b>	For use with Aquilion Short, Long, or One. Black, nylon, straps glide along slide rails on each side of table; affix to one another for easy securing of patients. Two at 6"	9061
9114	<b>Medium Bariatric Extension Strap</b>	For use with Aquilion Short, Long, or One. Affixes and lengthens 9062 Wide Security Strap for bariatric patients.	
9077	<b>Foot Extension Pad</b>	For use with Aquilion Short, Long, or One. Comfort Pad for Foot Extension, matches Table Pad construction and color, 1 x 18.33	9062
9112	<b>Foot Extension Pad Cover</b>	For use with Aquilion Short, Long, or One. Frosted vinyl pad cover, envelope-style. Flap tucks or slides under pad protecting joint, between couch and foot extension, from fluids.	
			9077




# Table Pads

Product Number	Product Description	Size	Application	
9055	Aquilion Short Table w/Protective Flaps	74.3" Long	Pad format makes it easier to move patients and adjust patient security straps which now slide along detached rail. Scan zone is indicated on protective flaps that keep fluid off of the table top and pallet.	 9056
9056	Aquilion Long Table w/Protective Flaps	86.1" Long	Pad format makes it easier to move patients and adjust patient security straps which now slide along detached rail. Scan zone is indicated on protective flaps that keep fluid off of the table top and pallet.	 9106
9071	Aquilion One Table w/Protective Flaps	96.75" Long	Pad format makes it easier to move patients and adjust patient security straps which now slide along detached rail. Scan zone is indicated on protective flaps that keep fluid off of the table top and pallet.	 9107
9106F	Aquilion Short Table w/Fluid Dam	73.8" Long		 9108
9107	Aquilion Long Table w/Fluid Dam	85.6" Long		
9108	Aquilion One Table w/Fluid Dam	96.3" Long		






# Table Pad Covers


Product Number	Product Description	Size	
9059	<b>Additional fluid protection of the scanner table and strapping system.</b> Secures to the topside of the fluid containment flaps with Velcoins.	Fits the Aquilion Short Table	
9060	<b>Additional fluid protection of the scanner table and strapping system.</b> Secures to the topside of the fluid containment flaps with Velcoins.	Fits the Aquilion Long Table	 9109
9072	<b>Additional fluid protection of the scanner table and strapping system.</b> Secures to the topside of the fluid containment flaps with Velcoins.	Fits the Aquilion One Table	
9109	<b>Frosted vinyl table cover "sleeve" with flap.</b> Access holes on bottom side of cover allow pad to adhere to table top Velcro. Compatible with all previous table pad designs.	Fits the Aquilion Short Table	 9110
9110	<b>Frosted vinyl table cover "sleeve" with flap.</b> Access holes on bottom side of cover allow pad to adhere to table top Velcro. Compatible with all previous table pad designs.	Fits the Aquilion Long Table	
9111	<b>Frosted vinyl table cover "sleeve" with flap.</b> Access holes on bottom side of cover allow pad to adhere to table top Velcro. Compatible with all previous table pad designs.	Fits the Aquilion One Table	 9111
<b>NEW!</b> 9121	<b>Frosted vinyl table cover fitting multiple table pads.</b>	Fits Aquilion Table Pads: 9106F, 9107 and 9108.	

# Detachable Side Rails

Product Number	Product Description	Size	
9078	<b>Velcro rails</b> attach to top and underside edges of Aquilion table top along which slide the 9061 and 9062 security straps.	Fits the Aquilion Short Table	 9078
9079	<b>Velcro rails</b> attach to top and underside edges of Aquilion table top along which slide the 9061 and 9062 security straps.	Fits the Aquilion Long Table	 9079
9080	<b>Velcro rails</b> attach to top and underside edges of Aquilion table top along which slide the 9061 and 9062 security straps.	Fits the Aquilion One Table	 9080

# Accessories

Product Number	Product Description	Size	Application	
9061	<b>Wide Security Strap</b>	14" Wide	Security straps slide along the detachable slide rail at the edge of the table pad to provide easy and comfortable patient security.	 9061
9062	<b>Narrow Security Strap</b>	6" Wide	Security straps slide along the detachable slide rail at the edge of the table pad to provide easy and comfortable patient security.	 9062
9113	<b>Wide Bariatric Extension Strap</b>	14" Wide	Affixes and lengthens part 9061 wide security strap for bariatric patients.	 9113
9114	<b>Narrow Bariatric Extension Strap</b>	6" Wide	Affixes and lengthens part 9062 narrow security strap for bariatric patients.	 9114
9001	<b>Rolled Edge Foot Extension Pad</b> Foam, Coated	1" x 14.17" x 18.33"	Comfort pad for table top foot extension. Fits all Aquilion series tables.	 9001

Product Number	Product Description	Size	Application	
9077	<b>Rolled Edge Foot Extension Pad</b> Compression Formed	1" x 14.17" x 18.33"	Comfort pad for table top foot extension. Fits all Aquilion series tables.	
9112	<b>Foot Extension Pad Cover</b>	Fits the 9077 Foot Extension Pad	Frosted vinyl pad cover, envelope style. Flap tucks or slides under pad protecting joint between couch and foot extension, from fluids.	 9112
9006	<b>Chin Strap</b>		Placed under the patient's chin to comfortably secure the patient in the head holder, allowing for stable and consistent imaging.	 9006
9007	<b>Forehead Strap with Adult Pad</b>		Placed around the patient's forehead to comfortably secure the patient in the head holder, allowing for stable and consistent imaging.	 9007
9008	<b>Adult Head Rest</b>	(1) Medium and (1) Large size	Comfort and stabilization pads for the CT Head Holder.	
9009	<b>Tilt Wedge</b>	9" x 6.5" x 3.38"	Unique, angled patient positioning wedge to assist in imaging procedures.	 9009
9017	<b>Knee Wedge</b>	21.88" x 18.33" 10.70"	Used to comfortably elevate the patient's knees to provide comfort for imaging procedures.	 9017
9025	<b>Coronal Head Extremity Holder</b>		This comfortable, prone-imaging, head holder provides unobstructed imaging of the temples. Also serves as a foot positioning device by placing the heel in the chin pad and the ball of the foot on the forehead pad.	 9025
9013	<b>Child Hugger with Table Adapter</b>	30.00" x 17.72" x 5.06"	Unique design safely and comfortably secures pediatric patients during CT scans.	
9014	<b>Forehead Strap with Child Pad</b>		Similar to adult head strap with child sized comfort positioning pad for use in the CT head holder.	 9013

Number	Product Description	Includes	Application
9015	<b>Child Head Cushion Kit</b>	(1) Small Size and (1) Medium Size	Comfortable positioning of child in CT head holder.
9019	<b>Supplemental Pediatric Accessory Kit</b>	(1) 9012 - Pediatric Hugger w/Table Adapter (1) 9013 - Child Hugger w/Table Adapter (1) 9014 - Forehead Strap w/Child Pad (1) 9015 - Child Head Cushion Kit	
9020F	<b>Pediatric Lift Pad</b> Compression Formed		Provides comfortable positioning for children whose shoulders rest on the attachment point of the head holder.












## Carbon Fiber Head Holder

Product Number	Product Description	Includes	Size
9090	<b>Aquilion Pediatric &amp; Adult Head Holder Kit</b>	(1) 9091 - Carbon Fiber Head Holder (1) 9092 - Adult Head Pad (1) 9093 - Child Head Pad (1) 9094 - Pediatric Head Pad (1) 9095 - Side Wedge Pads (set of 2) (1) 9096 - Bottom Wedge Pads (set of 2) (1) 9097 - Head Straps (set of 3)	Fits all Aquilion Systems
9091	<b>Carbon Fiber Head Holder</b>		
9092	<b>Adult Head Pad</b>		
9093	<b>Child Head Pad</b>		
9094	<b>Pediatric Head Pad</b>		
9095	<b>Side Wedges</b>	set of 2	
9096	<b>Bottom Wedges</b>	set of 2	
9097	<b>Head Straps</b>	set of 3	





# Veterinary Positioning Aids

Product Number	Product Description	Size	Application	
801-V	<b>Bolster</b>	8"H x 6.5"W x 16"L	Used for general head or body positioning for a wide variety of imaging views.	 801-V
802-V	<b>Body Cradle Large</b>	8"H x 16"W x 24"L	Radiolucent body cradles provide full-body positioning and support for both small and large breeds. Available with easily cleaned, 8-mil vinyl cover for protection against fluid contamination.	 802-V
802-V-8	<b>Body Cradle Large with 8-mil cover</b>	8"H x 16"W x 24"L	Radiolucent body cradles provide full-body positioning and support for both small and large breeds. Available with easily cleaned, 8-mil vinyl cover for protection against fluid contamination.	 802-V
803-V	<b>Body Cradle Small</b>	6"H x 12"W x 18"L	Radiolucent body cradles provide full-body positioning and support for both small and large breeds. Available with easily cleaned, 8-mil vinyl cover for protection against fluid contamination.	 803-V
803-V-8	<b>Body Cradle Small with 8-mil cover</b>	6"H x 12"W x 18"L	Radiolucent body cradles provide full-body positioning and support for both small and large breeds. Available with easily cleaned, 8-mil vinyl cover for protection against fluid contamination.	 803-V
804-V	<b>Skull Positioner Large</b>	6.5"H x 7.5"W x 10.5"L	Radiolucent skull positioners provide head positioning and support for small or large breeds for AP-PA lateral views.	 804-V
805-V	<b>Skull Positioner Small</b>	6"H x 6"W x 10"L	Radiolucent skull positioners provide head positioning and support for small or large breeds for AP-PA lateral views.	 804-V
812-V	<b>Rectangle</b>	1"H x 4"W x 8"L	Wedges and rectangles are used for corrective positioning of the skull, torso and extremities.	 816-V
816-V	<b>Bolster</b>	7"H x 10"W x 18"L	Used for general head or body positioning for a wide variety of imaging views.	 816-V



Product Number	Product Description	Size	Application
817-V	U-Block Pad	7.25"H x 14"W x 14"L	Used for general head or body positioning for a wide variety of imaging views.
818-V	Large U-Block	7.25"H x 14"W x 11"L	Used for general head or body positioning for a wide variety of imaging views.
822-V	Vet Echo Mattress	6"H x 18"W x 39"L	Mattress with cutout used for ultrasound procedures. Made with heavy, durable cover.
815-V	Triangle - 19° set of 2	3"H x 7"W x 8"L	Wedges and rectangles are used for corrective positioning of the skull, torso and extremities.
813-V	Rectangles set of 4	1"H x 4"W x 4"L	Wedges and rectangles are used for corrective positioning of the skull, torso and extremities.



## Veterinary Positioning Kits

Product Number	Product Description	Includes
808-V	General Vet Kit	(1) 117 - wedges (set of 2) (2) 118 - wedge (2) 139 - wedge (2) 126 - rectangles (set of 4) (1) 801-V
809-V	Small Animal Kit	(1) 105 - wedges (set of 2) (1) 113 - wedges (set of 2) (2) 126 - rectangles (set of 4)
810-V	Large Vet Set	(1) 117 - wedges (set of 2) (2) 118 - wedge (1) 122 - wedges (set of 2) (2) 139 - wedge (1) 801-V (1) 812-V (1) 813-V
814-V	Everyday Vet Set	(1) 116 - wedges (set of 2) (1) 154 - thin mattresses (set of 2) (1) 815-V (1) 816-V (1) 817-V (1) 818-V

# Urology Accessories

Product Number	Product Description	Units	Application	
21848	<b>Urology Hoop Frame and Brackets Unit</b> Fits U.S. & European Rails	EACH	Bracket and hoop frame integrated as one piece. Quick and easy installation.	 21848
27191	<b>Siemens Uroskop Table Pad</b>	EACH	Compression formed closed cell foam.	 27191

# Index

Part #	Description.....	Page	Part #	Description.....	Page	Part #	Description.....	Page
1	Deluxe Kit .....	11	126	Rectangle - Set of 4 .....	16	173	Table Pad - Foam Pad Only .....	19
2	Large Kit .....	11	127	Child AP-PA Head Holder .....	14	174	Table Pad - Foam Pad Only .....	19
3	Small Kit .....	11	128	Wedge - 13.5° .....	17	175	Table Pad - Foam Pad Only .....	19
4	Kit A.....	11	129	Wedge - 21°.....	17	176	Table Pad - Foam Pad Only .....	19
5	Kit B.....	11	130	Wedge - 45°.....	17	170-180	Foam Pad w/ 8 Standard Cover .....	19
6	Kit C.....	11	131	Rectangle .....	16	171-181	Foam Pad w/ 8 Standard Cover .....	19
7	Kit D.....	12	132	Head and Neck Support .....	9	172-182	Foam Pad w/ 8 Standard Cover .....	19
8	Kit E.....	12	133	Cervical Head Rest .....	9	173-183	Foam Pad w/ 8 Standard Cover .....	19
9	Rectangle Kit.....	12	135	EZ Straps - 1.5" x 30' roll - White hook & loop .....	20	174-184	Foam Pad w/ 8 Standard Cover .....	19
10	Kit F.....	12	136	Wedge - 12.5° .....	17	175-185	Foam Pad w/ 8 Standard Cover .....	19
11	Special Procedure Skull Kit .....	12	138	Wedge - 21°.....	17	176-186	Foam Pad w/ 8 Standard Cover .....	19
14	Ortho Kit .....	12	139	Wedge - 45°.....	18	191	Wedge - 45°.....	18
15	Kit I .....	13	140	Wedge - 27°.....	18	192	Leg Rest w/ Knee Cutouts – SM .....	21, 22
17	Routine X-Ray Kit .....	13	141	Spinal Wedge - 45° .....	18	193	Leg Rest w/ Knee Cutouts – M .....	21, 22
19	Standard R & F Kit .....	13	142	Rectangle .....	16	194	Leg Rest w/ Knee – L .....	21, 22
99	CT Kit .....	13, 22	143	Rectangle .....	16	195	Cardiac Wedge .....	7
100	Torso Block - 27°.....	6	144	Wedge - 45°.....	18	196	Wedge - 35/55°.....	18
101	Basal Block .....	9	145	Spinal Wedge - 45° .....	6, 18	199	Weight-Bearing Cassette Holder.....	20
102	AP-PA Head Holder.....	9	146	Cylinder.....	6	200	CT Hugger Complete Unit.....	14, 22
103	Rectangle .....	16	148	Child Myelography Block .....	15	201	Wedge - 35/55°.....	18
104	Pressure Block .....	9	149	Circular Disc .....	9	203	Cross-Table Lateral Leg Block .....	8
105	Wedge - 18.75° .....	16	153	Cylinder.....	6	204	Decubitus Block w/ Film Cassette Holder.....	20
106	Wedge - 45°.....	16	150	Myelography Block .....	6	205	Oblique Cervical Skull Block.....	9
107	Angular Bolster .....	6	151	Rectangle .....	16	207	4-in-1 Lateral Skull Cube.....	9
108	Wedge - 45°.....	16	152	Decubitus Block .....	16	208	Auto-Traction Straps .....	20
109	Angular Bolster .....	6	154	Thin Mats - Set of 2.....	16	209	Knee Arthrography Block .....	8
110	Wedge - 45° (set of 2) .....	17	156	CT Pillow .....	22	215	Endovaginal Wedge.....	20
111	Wedge - 45°.....	17	158	Wedge - 15° (set of 2) .....	18	250	Armboard Rail Mount .....	20
112	Circular Disc .....	9	159	Wedge - 20° (set of 2) .....	18	258	Armboard Pad - Fits 250 & 254.....	20
113	Wedge - 12° (set of 2) .....	17	160	Lumbar Support .....	6	254	Armboard Shoulder Mount .....	21
114	Rectangle .....	16	161	Lumbar Support .....	6	261	Rest Eazy Armboard GE .....	21
115	Bolster .....	6	162	Lumbar Support .....	6	262	Rest Eazy Armboard S&P .....	21
116	Wedge - 34° (set of 2) .....	17	163	Concave Disc .....	9	263	Rest Eazy Armboard Canon .....	21
117	Wedge - 16° (set of 2) .....	17	164	Concave Disc .....	9	265	Rest Eazy Armboard .....	21
118	Wedge - 45°.....	17	165	Concave Disc .....	9	300	Head Paddles - Packaged in Pairs .....	10
119	Head and Neck Rest.....	9	168	Security Strap .....	20	301	Pediatric Dual Oblique Finger Block .....	15
120	Wedge - 18°.....	17	169	Decubitus Block .....	16	302	Pediatric Dual Lateral Finger Block .....	15
121	Oblique Finger Block.....	7	170	Table Pad - Foam Pad Only .....	19	303	Head-Ends - Packaged in Pairs .....	10
122	Wedge - 22° (set of 2) .....	17	171	Table Pad - Foam Pad Only .....	19	304	AP-PA Head Immobilizer - Pediatric .....	15
123	Wedge - 33°.....	17	172	Table Pad - Foam Pad Only .....	19	305	AP-PA Head Immobilizer - Child.....	15
124	Bolster .....	6				308	Pediatric AP-PA Hand/Foot I	

# Index

Part #	Description	Page	Part #	Description	Page	Part #	Description	Page
	Immobilizer	15	9001	Rolled Edge Foot Extension Pad – Foam, Coated	30	9094	Pediatric Head Pad	32
401	Leg Immobilizer	8	9006	Chin Strap	31	9095	Side Wedges - set of 2	32
405	Waters Positioner	10	9007	Forehead Strap with Adult Pad	31	9096	Bottom Wedges - set of 2	32
407	Dual Lateral Skull Block	10	9008	Adult Head Rest	31	9097	Head Straps - set of 3	32
408	AP-PA Head Immobilizer	10	9009	Tilt Wedge	31	9100	Aquilion Standard Kit - Short	23
409	AP Head Immobilizer	10	9013	Child Hugger w/ Table Adapter	31	9101	Aquilion Standard Kit - Long	24
422	Block Hanger Kit	14	9014	Forehead Strap w/ Child Pad	31	9102	Aquilion Standard Kit - Aquilion One	25
423	Velcoins - set of 24	14	9015	Child Head Cushion Kit	32	9105F	Identical to 9102 above, but without the 9090 Head Holder Kit	25
425	The Pediatric Hugger Complete Unit	14	9017	Knee Wedge	31	9106F	Aquilion Short w/ Fluid Dam	28
426	The Child Hugger Complete Unit	14	9019	Supplemental Pediatric Accessory Kit	32	9107	Aquilion Long w/ Fluid Dam	28
427	Leg Brace	8	9020F	Pediatric Lift Pad - Compression Formed	32	9108	Aquilion One w/ Fluid Dam	28
428	Lateral Finger Block - Right Hand	7	9025	Coronal Head Holder	10, 31	9109	Frosted vinyl table cover “sleeve” w/ flap	29
429	Lateral Finger Block - Left Hand	7	9055	Aquilion Short w/ Protective Flaps	28	9110	Frosted vinyl table cover “sleeve” w/ flap	29
430	Dual Oblique Finger Block	7	9056	Aquilion Long w/ Protective Flaps	28	9111	Frosted vinyl table cover “sleeve” w/ flap	29
431	Dual Lateral Finger Block	7	9059	Additional fluid protection of the scanner table and strapping system	29	9112	Foot Extension Pad Cover	31
438	Rounded Shoulder Support	7	9060	Additional fluid protection of the scanner table and strapping system	29	9113	Wide Bariatric Extension Strap	30
439	Winged Shoulder Support	7	9061	Wide Security Strap	30	9114	Narrow Bariatric Extension Strap	30
463	MRI Comfort Pillow	21	9062	Narrow Security Strap	30	9121	Frosted vinyl table cover fitting multiple table pads	29
466	Sloped Leg Immobilizer	8	9071	Aquilion One w/ Protective Flaps	28	21848	Urology Hoop Frame & Brackets Unit - Fits U.S. & European Rails	35
520	Table Pad w/ heavy duty 20 mil vinyl cover	19	9072	Additional fluid protection of the scanner table and strapping system	29	27191	Siemens Uroskop Table Pad	35
530	Table Pad w/ heavy duty 20 mil vinyl cover	19	9077	Rolled Edge Foot Extension Pad - Compression Formed	31	27625	Manual Pump	22
801-V	Bolster	33	9078	Velcro rails attach to top and underside edges of Aquilion table top along which slide the 9061 and 9062 security straps	30	VMR33PUBB	Pediatric Immobilizer	14, 22
802-V	Body Cradle - Large	33	9079	Velcro rails attach to top and underside edges of Aquilion table top along which slide the 9061 and 9062 security straps	30	VMR38X01	Child Immobilizer	14, 22
802-V-8	Body Cradle - Large w/ 8-mil cover	33	9080	Velcro rails attach to top and underside edges of Aquilion table top along which slide the 9061 and 9062 security straps	30			
803-V	Body Cradle - Small	33	9090	Aquilion Pediatric & Adult Head Holder Kit	32			
803-V-8	Body Cradle - Small w/ 8-mil cover	33	9091	Carbon Fiber Head Holder	32			
804-V	Skull Positioner - Large	33	9092	Adult Head Pad	32			
805-V	Skull Positioner - Small	33	9093	Child Head Pad	32			
808-V	General Vet Kit	34						
809-V	Small Animal Kit	34						
810-V	Large Vet Set	34						
812-V	Rectangle	33						
814-V	Everyday Vet Set	34						
816-V	Bolster	33						
817-V	U-Block Pad	34						
818-V	Large U-Block	34						
822-V	Vet Echo Mattress	34						
815-V	Triangle - 19° - set of 2	34						
813-V	Rectangles - set of 4	34						



Distributed by



(877) 338-6857  
info@ultraray.com  
ultraraymedical.com