



**VETERINARY:**

## Radiation Protection and Accessories

Learn more at [ultraraymedical.com](http://ultraraymedical.com)

• RADIOLOGY

• CT and MRI

• ELECTROMEDICAL

• MAMMOGRAPHY

• NUCLEAR MEDICINE

• ULTRASOUND

• ONCOLOGY/RADIATION THERAPY

• STERILE PROCESSING

• OR SUPPLIES

• PATIENT HANDLING

• FURNITURE



Standard Thyroid Shield  
(Loose, Tethered or Buckle)



Integrated Thyroid Shield  
(Sewn-On)

## VA1 ADJUSTABLE FIT FRONTAL with Velcro Closure

- Easy on and off convenience.
- Adjustable velcro arms provide excellent weight distribution and comfort.
- Built-in shoulder pads offer additional comfort.
- Adjustable shoulder strap prevents the apron from slipping.

## VA2 ADJUSTABLE FIT FRONTAL with Belt and Buckle

- Easy on and off convenience.
- Adjustable heavy-duty 2" web belt with quick-release buckle closure.
- Built-in shoulder pads offer additional comfort.
- Adjustable shoulder strap prevents the apron from slipping.

## VA3 ECONOMY FRONTAL with Tie-Strap Closure

- Easy on and off convenience.
- 1" straps cross the back and tie securely at the front.

CODE (APRON)	DESCRIPTION
VA1	Veterinary Frontal with Velcro
VA2	Veterinary Frontal with Belt & Buckle
VA3	Veterinary Frontal with Tie-Straps

SIZE	ACTUAL APRON WIDTH X LENGTH	SHOULDER WIDTH	MAX. CHEST/HIP
X-Small	20" x 34"	14"	34"
Small	22" x 35"	14 <sup>1/2</sup> "	38"
Medium	24" x 37"	16"	42"
Large	24" x 39"	18 <sup>1/2</sup> "	44"
X-Large	26" x 41"	20 <sup>1/2</sup> "	46"
XX-Large	28" x 43"	21 <sup>1/2</sup> "	50"

CODE (THYROID SHIELD)	DESCRIPTION
TSR-IT	Thyroid Shield - Sewn-on (Integrated)
TSRL	Standard Thyroid Shield Loose
TSRT	Standard Thyroid Shield Tethered
TSRB	Standard Thyroid Shield Tethered with Buckles

## UltraLite

- Stylish and comfortable wraparound frame.
- High-impact resistant plastic for safety.
- Colours: Black, Bronze, Carbon Fibre, Pink, Purple/Silver, Red, Silver and Yellow/Blue.



## Sport Wrap

- Solid, sturdy and extremely comfortable wraparound glasses.
- Built-in side-shields for safety.
- Colours: Black/Yellow and Brown/Orange.
- Rx: Single, Bifocal and Progressive.



## Maxi (Mini, Regular or Large)

- Padded temple bars and nose bridge for extra comfort.
- The most durable product in our line.
- Colours: Black, Blue, Red and Silver.
- Rx: Single, Bifocal and Progressive.



## Basic 52 Wrap

- Excellent peripheral vision with leaded glass side-shields.
- Economical and lightweight with a spring hinge.
- Silicone, adjustable non-slip uni-fit bridge for increased comfort.
- Colours: Tortoise.
- Rx: Single, Bifocal and Progressive.



## Fitover (Regular or Large)

- Designed to fitover existing prescription glasses.
- Black rubberized nylon frame.
- Lightweight and comfortable.
- Colours: Black and Tortoise.

Regular



Large



\*All glasses are 0.75mm L.E. Side-Shields (If equipped) are 0.35mm L.E. All glasses include a protective case and strap.

## VET Mittens

This shield offers lead-equivalency of 0.5mm while positioning animals during procedures. Combination of Velcro wrist strap and non-slip material on palm prevents slippage while securing your patient.



## Open Hand Shield

This shield resembles a paddle with a cuff. It is used to protect the top of the hand. One-size fits all with a lead-equivalency of 0.5mm.



## Open Palm Mitts

These mitts are primarily designed for veterinary applications where fingers and thumbs are available to be used in a procedure. One-size fits all with a lead-equivalency of 0.5mm.



## Open Palm with Thumb

Identical to the Open-Palm Mitt with a closed thumb. One-size fits all with a lead-equivalency of 0.5mm.



## Ducky Mitts

These mitts are designed for veterinary applications where fingers and thumbs must be protected. One-size fits all with a lead-equivalency of 0.5mm.



## Closed Mitts with Thumb

Similar to the Ducky Mitt without the finger pouch. One-size fits all with a lead-equivalency of 0.5mm.



## 15" Lead Vinyl Gloves

Each glove features an outer one-piece soft urethane molded design with an inner liner made from breathable nylon for comfort and is easily removed for inspection. Available in pairs with a lead-equivalency of 0.5mm.





## VETERINARY SANDBAGS with/without handles

Sandbags are rectangular in shape. An inner bladder is filled with sand and then inserted into a blue nylon or vinyl bag.

CODE (NYLON)	CODE (VINYL)	DESCRIPTION
VSB5-N	VSB5-V	VET Sandbag - 5lb - 7" x 8"
VSB7-N	VSB7-V	VET Sandbag - 7lb - 7" x 12"
VSB12-N	VSB12-V	VET Sandbag - 12lb - 7" x 16"
VSB5H-N	VSB5H-V	VET Sandbag - 5lb - 7" x 8" with Handles
VSB7H-N	VSB7H-V	VET Sandbag - 7lb - 7" x 12" with Handles
VSB12H-N	VSB12H-V	VET Sandbag - 12lb - 7" x 16" with Handles
VSB12S-N	VSB12S-V	VET Sandbag - Snake - 12lb - 5" x 27"



## CUSTOM EMBROIDERY

- Custom embroidery and stock designs available.
- A logo or design specific to your facility will further enhance the product and clearly identify it.
- For more options, please refer to our main Medical Catalogue.



CANADA'S TRUSTED  
SOURCE FOR IMAGING  
PRODUCTS

