

TechnoAide

Distributed by



(877) 338-6857
info@ultraray.com
ultraraymedical.com

PROTECTIVE APPAREL & SUPPLIES

7

SPONGES & POSITIONING AIDS

20

TABLE PADS

35

X-RAY MARKERS

36

GENERAL SUPPLIES & ACCESSORIES

38

CLOSED CELL SPONGES

Techno-Aide's NEW closed cell sponges are the ultimate option in positioning aids! The liquid resistant and artifact free material offers every advantage of traditional coated sponges with ZERO of the traditional drawbacks! Durable, lightweight and cleanable, these non-porous positioning aids offer completely artifact free imaging without the need for expensive coating! Techno-Aide's closed cell sponges are the future of positioning aids.



**CLEANABLE • LIQUID PROOF
ARTIFACT FREE**

VINYL COVERED SPONGES

Techno-Aide's premium vinyl covered sponges offer superior durability and stability during extended imaging. These extremely durable and easy to clean vinyl covered sponges feature soft and flexible foam interiors for increased patient comfort. Techno-Aide's unique positioners feature completely sewn and sealed seams and edges, creating a totally sanitary and non-porous surface to prevent liquid absorption and cross-contaminating between patients.



COATED & NON-COATED SPONGES

Techno-Aide's coated and non-coated sponges offer excellent positioning support at an affordable price. High quality and durable foam provides added patient comfort and helps achieve perfect positioning during imaging. All Techno-Aide coated sponges feature a sprayed-on proprietary non-porous, cleanable and flexible green exterior.



Visit [TECHNO-AIDE.com](https://www.techno-aide.com) to shop our entire selection of Sponges & Positioning Aids!

A BRIEF GUIDE TO TECHNO-AIDE POSITIONING SPONGES

Non-Coated




Example - YFBT

\$

- Basic (core) piece of porous (open cell) foam
- No cleanable surface
- Not resistant to fluids
- Not considered durable
- Radiolucent
- Typically artifact free
- Most inexpensive option on average
- Some shapes available in "Stealth" design to reduce the chance of artifacts

Coated



Example - YCBT

\$\$\$

- Open cell foam coated with proprietary medical-grade coating
- Cleanable surface
- Resistant to fluids
- More durable than basic (core) foam
- Radiolucent
- Typically not artifact free
- Second most expensive option on average
- Some shapes available in "Stealth" design to reduce the chance of artifacts

Vinyl Covered



Example - YVBT

\$\$\$\$

- Open cell foam covered with completely sewn heavy-duty vinyl cover
- Cleanable surface
- Resistant to fluids
- Extreme durability and comfort
- Radiolucent
- Not artifact free (sewn seams will show during most imaging)
- Most expensive option on average
- Not available in "Stealth" design

Closed Cell



Example - YXBT

\$\$

- High density, crosslink (closed cell) foam
- Cleanable surface
- Entirely resistant to fluids
- Extreme durability
- Radiolucent
- 100% artifact free
- No coating/cover required
- Second most inexpensive option on average
- "Stealth" design not necessary

UNDERSTANDING "STEALTH" SPONGES

Non-Stealth or "Traditional" Sponges



Sponges formed and cut into traditional shapes, typically featuring standard right angles. This type is the most common style across all positioning sponges. Right angles typically show as artifacts during imaging, especially when a form of coating or covering is applied to the core sponge.

Stealth or "Flared" Sponges



Traditional sponge shapes reformed with proprietary "flared" sides, removing most right angles. This patented and proprietary design helps significantly reduce the chance of artifacts during imaging. Flared sides increase the overall size of the sponge, typically at the base of the product.

SPONGES & POSITIONING AIDS

CLOSED CELL SPONGES

**Measurements in L x W x H

FLAT BOTTOM SPONGES

YXAM - Flat Bottom V-Trough
Sponge [15" x 14" x 4.625"]

YXCM - LG Flat Bottom
V-Trough/AP & PA Adult Head
Holder Sponge [36" x 14" x 4.625"]



YXAM

EXTREMITY SPONGES

YXCA - S Extremity Sponge
[13" x 7" x 3.5"]

YXCB - M Extremity Sponge
[19.5" x 11.5" x 5"]

YXCC - L Extremity Sponge
[30.5" x 11" x 7"]

YXCD - XL Extremity Sponge
[44" x 14" x 9"]



YXCA

RECTANGLE SPONGES

YXBL - Rectangle Sponge
[10" x 7.5" x 3"]

YXBN - Rectangle Sponge
[18" x 14" x 3"]

YXBR - Rectangle Sponge
[6" x 4" x 1"]

YXBT - Rectangle Sponge
[6" x 4" x 2"]

YXBV - Square Block Sponge
[6" x 6" x 2"]

YXGF - Rectangle Sponge
[12" x 10" x 2"]

YXMB - Rectangle Sponge
[10" x 8" x 2"]



YXBL

WEDGE SPONGES

YXBE - 34° Wedge Sponge
[7" x 4" x 3"]

YXBF - 27° Wedge Sponge
[10" x 10" x 6"]

YXBG - 21° Wedge Sponge
[10" x 10" x 4.625"]

YXBI - 15° Wedge Sponge
[10" x 7" x 2"]

YXBJ - 12° Wedge Sponge
[9.5" x 7.25" x 2.25"]



YXBE

FINGER BLOCK SPONGES

YXAA - Finger Block Sponge
[5" x 3" x 3"]

YXFL - L Lateral Finger Block
Sponge [6" x 6.125" x 3.25"]

YXFR - R Lateral Finger Block
Sponge [6" x 6.125" x 3.25"]

YXAD - Adult Bi-Lateral Finger
Block Set [YXFL + YXFR]



YXAA

45° WEDGE SPONGES

YXBB - M 45° Spinal Body
Wedge Sponge [14" x 7" x 7"]

YXBC - 45° Wedge Sponge
[7" x 7" x 7"]

YXBD - LG 45° Spinal Body
Wedge Sponge [21" x 7" x 7"]

YX16 - 45° Wedge Sponge
[30" x 9" x 9"]



YXBB

CLOSED CELL SPONGES

**Measurements in L x W x H

SPECIALTY SPONGES



YXBU – 30° – 60° – 90° Multiangle Wedge [29.5" x 10" x 6"]



YXAY – Specialty Angular Bolster Sponge [14" x 22.5" x 8"]



YXWB – Weight Bearing Cassette Holder [17.5" x 11.5" x 4" | Channel: 1.1875" x 11.5" x 3"]



YXAJ – Knee Wedge Sponge [22" x 13.5" x 7"]



YXAK – Circular Head Sponge [Diameter: 7" x 1.5"]



YXNZ – Podiatry Axial/Sesamoid Weight Bearing Sponge [13" x 8" x 4"]



YXUZ – Ultrasound Wedge Sponge [20" x 20" x 6"]

CLOSED CELL BUNDLES


| | |
|----------|--|
| YZCA | Clinic Sponge Bundle A (1 each – YXAA, YXAK, YXBB, YXBI) |
| YZCB | Clinic Sponge Bundle B (1 each – YXAA, YXAK, YXBC, YXBD, YXBE, YXBI, YXBJ, YXBL) |
| YZCC | General Sponge Bundle 1 (1 each – YXAA, YXAK, YXBD, YXBI, YXBN, YXBV, YXMB) |
| YZGA | General Sponge Bundle A (1 each – YXAK, YXBC, YXBE, YXBI, YXBL) |
| YZGB | General Sponge Bundle B (1 each – YXAA, YXAK, YXBC, YXBE, YXBI, YXBJ, YXBL) |
| YZGP | General Sponge Bundle P (1 each – YXAA, YXAK, YXBD, YXBI) |
| YB-X-A01 | Veterinary Sponge Bundle (1 each – YXBI, YXBT, YXBR, YXAM, YXBD) |
| YB-X-A02 | Urgent Care Sponge Bundle (1 each – YXAA, YXAD, YXAK, YXBC, YXBD, YXBE, YFBJ, YXBL, YXBN) |
| YB-X-A03 | Alt. Urgent Care Sponge Bundle (1 each – YXAA, YXAD, YXAK, YXBC, YXBD, YXBE, YXBI, YXBL, YXBN) |

SPONGES & POSITIONING AIDS

COATED SPONGES

**Measurements in L x W x H


45° SPONGES

| | | |
|---|---|--|
| YCBA – Bariatric Spinal Wedge Sponge (Stealth) [27.5" x 10" x 7"] |  <p>YCBA</p> | |
| YCBB – Spinal Wedge Sponge (Stealth) [20.25" x 10" x 7"] | | |
| YCBC – Spinal Body Wedge Sponge (Stealth) [13.5" x 10" x 7"] | | |
| YCBD – Spinal Body Wedge Sponge (Stealth) [27.75" x 10" x 7"] | | |
| YCCG – Bariatric Body Wedge Sponge (Non-Stealth) [30" x 9" x 9"] | | YCEB – Wedge Sponge (Stealth) [30" x 10" x 7.5"] |
| YCDA – Cardiac Wedge Sponge (Non-Stealth) [28" x 19" x 14"] | | YCEE – Echocardiography Wedge Sponge (Stealth) [38.25" x 13" x 9"] |
| YCEA – Spinal Wedge Sponge (Stealth) [33" x 13.75" x 9.5"] | YCFD – Spinal Body Wedge Sponge (Non-Stealth) [21" x 7" x 7"] | |


WEDGE SPONGES

| | |
|---|---|
| YCAJ – 26° Knee Wedge Support Sponge (Stealth) [21" x 19" x 8"] |  <p>YCBF</p> |
| YCBE – 34° Wedge Sponge (Stealth) [9.75" x 5.5" x 3"] | |
| YCBF – 27° Wedge Sponge (Stealth) [15.5" x 12.75" x 6"] | |
| YCBG – 21° Wedge Sponge (Stealth) [14.25" x 12.25" x 4.5"] | YCBU – 30°-60°-90° Multiangle Wedge (Stealth) [29.5" x 10" x 6"] |
| YCBH – 16° Wedge Sponge Set (Stealth) [11" x 8.5" x 3.5"] | YCDB – 30° Wedge Sponge (Non-Stealth) [27" x 28" x 16"] |
| YCBI – 15° Wedge Sponge Set (Stealth) [10" x 12.5" x 3"] | YCEG – 65° Wedge Sponge (Non-Stealth) [25" x 16" x 8.25"] |
| YCBJ – 12° Wedge Sponge Set (Stealth) [9.25" x 10.5" x 2.25"] | YCEU – 30°-60°-90° Multiangle Bariatric Wedge Sponge (Stealth) [29.5" x 10" x 6"] |

RECTANGLE SPONGES

| | |
|--|---|
| YCBL – Rectangle Sponge (Stealth) [12.5" x 10.25" x 3"] |  <p>YCBL</p> |
| YCBM – Bariatric Rectangle Sponge (Stealth) [14.5" x 10.5" x 3"] | |
| YCBN – Rectangle Sponge (Stealth) [20.5" x 16.25" x 3"] | |
| YCBO – Square Sponge (Stealth) [13.5" x 13.5" x 4"] | YCBT – Rectangle Sponge (Stealth) [6" x 4" x 2"] |
| YCBP – Bariatric Decubitus Rectangle Sponge (Stealth) [27.5" x 21.5" x 4"] | YCBV – Square Sponge (Non-Stealth) [6" x 6" x 2"] |
| YCBQ – Rectangle Sponge (Non-Stealth) [6" x 4" x .5"] | YCBY – Rectangle Sponge (Stealth) [27.5" x 21.5" x 4"] |
| YCBR – Rectangle Sponge (Non-Stealth) [6" x 4" x 1"] | YCEM – Rectangle Sponge (Non-Stealth) [24" x 15" x 2.5"] |
| YCBS – Rectangle Sponge (Non-Stealth) [6" x 4" x 1.5"] | YCGG – Rectangle Sponge (Stealth) [19.25" x 9.25" x 1.5"] |

HEAD SPONGES

| | |
|---|---|
| YCAK – Circular Head Block (Stealth) [Diameter: 7" (top) x 19" (base) x 1.5"] |  <p>YCAK</p> |
| YCAL – Head & Neck Support Sponge (Stealth) [14" x 16" x 6.5"] | |
| YCAM – AP/PA Adult Head Holder Sponge (Stealth) [15" x 14" x 6.75"] | |
| YCAN – AP/PA Child Head Holder Sponge (Stealth) [11" x 9.75" x 5"] | YCBW – CT Head Rest Sponge (Non-Stealth) [12.5" x 11" x 10"] |
| YCAP – Adult Head Immobilizer Sponge (Stealth) [15.25" x 9" x 7"] | Y CET – Occipital Head Support Sponge (Non-Stealth) [13" x 12.75" x 4.5"] |
| YCAQ – Child Head Immobilizer Sponge (Stealth) [11.5" x 6.75" x 5.25"] | YCEV – CT/PET Head Holder Sponge (Non-Stealth) [13.5" x 14" x 6.75"] |
| YCAR – Cervical Head Rest Sponge (Stealth) [16" x 7.5" x 4.5"] | YCEW – Philips CT Head Rest Sponge (Non-Stealth) [12" x 12.5" x 6.5"] |

COATED SPONGES

**Measurements in L x W x H

EXTREMITY SPONGES

| |
|--|
| YCCA - SM Extremity Sponge (Stealth) [12.5" x 7" x 3"] |
| YCCB - MD Extremity Sponge (Stealth) [19.75" x 11.5" x 5"] |
| YCCC - LG Extremity Sponge (Stealth) [30.5" x 11" x 7"] |
| YCCD - XL Extremity Sponge (Stealth) [44.5" x 14" x 9"] |



YCCA

ULTRASOUND SPONGES

| |
|---|
| YCEP - Bariatric Endo-Ultrasound Wedge Sponge (Non-Stealth) [16" x 16" x 4.5"] |
| YCEQ - Perineal Endovaginal/Transvaginal Ultrasound Wedge Sponge (N-S) [20" x 20" x 3"] |
| YCUZ - Endo-Ultrasound Wedge Sponge (Non-Stealth) [20" x 20" x 6"] |



YCEP

FINGER SPONGES

| | |
|---|---|
| YCAA - Adult Oblique Finger Block Sponge (Stealth) [10.5" x 6" x 3"] | |
| YCAB - Bi-lateral Adult Oblique Finger Block Sponge (Stealth) [12.5" x 10.5" x 3"] | |
| YCAC - Bi-lateral Pediatric Oblique Finger Block Sponge (Stealth) [7.75" x 7.5" x 2"] | YCAA |
| YCAD - Adult Lateral Finger Block Sponge (Stealth) [13" x 8" x 3"] | YCFR - R Lateral Finger Block Sponge [5-6.5" x 7.5" x 3"] |
| YCAE - Pediatric Lateral Finger Block Sponge [7" x 4" x 1.75"] | YCFL - L Lateral Finger Block Sponge [5-6.5" x 7.5" x 3"] |

LEG SPONGES

| |
|---|
| YCAG - Bi-lateral Leg Immobilizer Sponge (Stealth) [21" x 20" x 6"] |
| YCAH - Knee Anthrogram Block Sponge (Stealth) [19" x 14.5" x 8"] |
| YMCM - MD Leg Rest w/ Knee Cutouts (Stealth) [16.25" x 12.5" x 7"] |
| YMCP - LG Leg Rest w/ Knee Cutouts (Stealth) [24" x 18" x 10"] |



YCAG

YMCS - SM Leg Rest w/ Knee Cutouts (Stealth) [18" x 14.5" x 9"]


SPECIALTY SPONGES



YCAT - Bariatric Torso Sponge (Stealth) [16.5" x 18.5" x 7"]



YCAW - Adult Myelogram Block Sponge (Non-Stealth) [22" x 10" x 14"]



YCFQ - Decubitus Block w/ Imaging Plate Holder (Non-Stealth) [Lower: 36" x 16" x 4" | Back: 36" x 10" x 11" | Channel: 36" x .05" x 3"]



YCAV - Specialty Angular Bolster Sponge (Stealth) [22.5" x 14" x 7.25"]



YCAX - Child Myelogram Block Sponge (Non-Stealth) [14" x 8" x 10"]



YCMG - Contoured Torso Sponge (Non-Stealth) [18.5" x 17" x 5"]

SPONGES & POSITIONING AIDS

COATED SPONGES

**Measurements in L x W x H

COATED BUNDLES

| | |
|----------|--|
| YSAS | Knee Arthropasty Stress Assembly Sponge Set (Non-Stealth) |
| YSBA | Bariatric Sponge Bundle (1 each - YCAT, YCBA, YCBM, YCBP) |
| YSCA | Clinic Sponge Bundle A (1 each - YCAA, YCAK, YCBB, YCBI) |
| YSCB | Clinic Sponge Bundle B (1 each - YCAA, YCAK, YCBC, YCBD, YCBE, YCBI, YCBJ, YCBL) |
| YSCT | CT Sponge Bundle A (1 each - YCBB, YCBC, YCBE, YCBF, YCBG, YCBI, YCBJ, YCBW) *Bundle non-returnable. More info upon request. |
| YSGA | General Sponge Bundle A (1 each - YCAK, YCBC, YCBE, YCBJ, YCBL) |
| YSGB | General Sponge Bundle B (1 each - YCAA, YCAK, YCBC, YCBE, YCBI, YCBJ, YCBL) |
| YSGC | General Sponge Bundle C (1 each - YCAA, YCAM, YCAR, YCAY, YCBL, YCAH) |
| YSGD | General Sponge Bundle D (1 each - YCAL, YCAN, YCAR, YCBF, YCBJ, YCBL) |
| YSGE | General Sponge Bundle E (1 each - YCAA, YCAK, YCAM, YCBB, YCBG, YCBI, YCBL, YCBN, YCBU) |
| YSGG | General Sponge Bundle G (1 each - YCAA, YCAK, YCBD, YCBE, YCBJ, YCBL) |
| YSGP | General Sponge Bundle P (1 each - YCAA, YCAK, YCBD, YCBI) |
| YSOA | Orthopedic Sponge Bundle A (1 each - YCAA, YCAM, YCAR, YCBD, YCBL) |
| YSOB | Orthopedic Sponge Bundle B (1 each - YCAL, YCAN, YCAR, YCBC, YCBG, YCBJ, YCBL) |
| YSOD | Orthopedic Sponge Bundle D (1 each - YCAA, YCAR, YCBD, YCBI, YCBL) |
| YSHA | Head & Neck Sponge Bundle (1 each - YCAK, YCAL, YCAM, YCAN, YCAP, YCAQ, YCAR) |
| YSPA | Pediatric Sponge Bundle (1 each - YCAC, YCAE, YCAK, YCAN, YCAQ, YCAX) |
| YSRA | Rectangle Sponge Bundle (3 each - YCBQ, YCBR, YCBS, YCBT) |
| YSUB | Ultrasound Breast or Paracentesis Sponge Bundle A (1 each - YCBG, YCFD, YCBU) |
| YSVA | General Sponge Bundle 1 (2 each - YCAK; 1 each - YCBD, YCBI, YCBL, YCBT) |
| YSWA | Wedge Sponge Bundle A (1 each - YCBC, YCBF, YCBG, YCBI) |
| YSWB | Wedge Sponge Bundle B (1 each - YCBD, YCBE, YCBF, YCBG) |
| YB-C-A02 | Urgent Care Sponge Bundle (1 each - YCAA, YCAD, YCAK, YCBC, YCED, YCBE, YCBJ, YCBL, YCBN) |
| YSRC | Scanning Support Sponge Bundle (10 each - YCBT) |
| YSRD | Ultrasound Support Sponge Bundle (1 each - YCAY, YCEE, YCUZ) |


NON-COATED SPONGES

**Measurements in L x W x H


45° WEDGE SPONGES

| | |
|--|---|
| YBBA - Bariatric Spinal Wedge Sponge (Stealth) [27.5" x 10" x 7"] |  <p>YFBB</p> |
| YFBB - Spinal Wedge Sponge (Stealth) [20.25" x 10" x 7"] | |
| YFBC - Spinal Body Wedge Sponge (Stealth) [13.5" x 10" x 7"] | |
| YFBD - Spinal Body Wedge Sponge (Stealth) [27.75" x 10" x 7"] | |
| YFEA - Spinal Wedge Sponge (Stealth) [24" x 13.75" x 9.5"] | |
| YFEE - Echocardiography Wedge Sponge (Stealth) [38.25" x 13" x 9"] | YFFD - Spinal Wedge Sponge (Non-Stealth) [21" x 7" x 7"] |

WEDGE SPONGES

| | |
|---|---|
| YFBE - 34° Wedge Sponge Set (Stealth) [9.75" x 5.5" x 3"] |  <p>YFBE</p> |
| YFBF - 27° Wedge Sponge (Stealth) [15.5" x 12.75" x 6"] | |
| YFBG - 21° Wedge Sponge (Stealth) [14.25" x 12.25" x 4.5"] | |
| YFBH - 16° Wedge Sponge Set (Stealth) [11" x 8.5" x 3.5"] | |
| YFBI - 15° Wedge Sponge Set (Stealth) [10" x 12.5" x 3"] | |
| YFBJ - 12° Wedge Sponge Set (Stealth) [9.25" x 10.5" x 2.25"] | YFBU - 30°-60°-90° Multiangle Wedge Sponge Set (Stealth) [29.5" x 10" x 6"] |
| YFDB - 30° Cardiac Wedge Sponge Set (Stealth) [27" x 28" x 16"] | |


RECTANGLE SPONGES

| | |
|---|--|
| YFBL - Rectangle Sponge (Stealth) [12.5" x 10.25" x 3"] |  <p>YFBO</p> |
| YFBN - Rectangle Sponge (Stealth) [20.5" x 16.25" x 3"] | |
| YFBO - Square Sponge (Stealth) [13.5" x 13.5" x 4"] | |
| YFBR - Rectangle Sponge (NonStealth) [6" x 4" x 1"] | |
| YFBT - Rectangle Sponge (Stealth) [6" x 4" x 2"] | |
| YFBY - Rectangle Sponge (Stealth) [27.5" x 21.5" x 4"] | YFEM - Rectangle Sponge (Non-Stealth) [24" x 15" x 2.5"] |


HEAD SPONGES

| | |
|--|---|
| YFAK - Circular Head Block (Stealth) [Diameter: 7" (top) - 9" (base) - 1.5"] |  <p>YFAK</p> |
| YFAL - Head & Neck Support Sponge (Stealth) [16" x 14.5" x 6.5"] | |
| YFAM - AP/PA Adult Head Holder Sponge (Stealth) [15" x 14" x 6.75"] | |
| YFAN - AP/PA Child Head Holder Sponge (Stealth) [11" x 9.75" x 5"] | |
| YFAQ - Child Head Immobilizer Sponge (Stealth) [11.5" x 6.75" x 5.25"] | |
| YFAR - Cervical Head Rest Sponge (Stealth) [16" x 7.5" x 4.5"] | YFBW - CT Head Rest Sponge (Non-Stealth) [12.5" x 11" x 10"] |

EXTREMITY SPONGES

| | |
|--|---|
| YFCA - SM Extremity Sponge (Stealth) [12.5" x 7" x 3"] |  <p>YFCA</p> |
| YFCB - MD Extremity Sponge (Stealth) [19.75" x 11.5" x 5"] | |
| YFCC - LG Extremity Sponge (Stealth) [24" x 11" x 7"] | |
| YFCD - XL Extremity Sponge (Stealth) [34" x 14" x 9"] | |

FINGER SPONGES

| | |
|--|---|
| YFAA - Adult Oblique Finger Block Sponge (Stealth) [10.5" x 6" x 3"] |  <p>YFAB</p> |
| YFAB - Bi-lateral Adult Oblique Finger Block Sponge (Stealth) [12.5" x 10.5" x 3"] | |
| YFAC - Adult Lateral Finger Block Sponge (Stealth) [7.75" x 7.5" x 2"] | |
| YFAD - Adult Lateral Finger Block Sponge (Stealth) [13" x 8" x 3"] | |
| YFAE - Pediatric Lateral Finger Block Sponge [7" x 4" x 1.75"] | |

SPONGES & POSITIONING AIDS

NON-COATED SPONGES

**Measurements in L x W x H

SPECIALTY SPONGES



YFAW - Adult Myelogram Block Sponge (Non-Stealth) [22" x 10" x 14"]



YFAX - Child Myelogram Block Sponge (Non-Stealth) [14" x 8" x 10"]



YFAY - Specialty Angular Bolster Sponge (Stealth) [22.5" x 14" x 7.25"]

NON-COATED BUNDLES

| | |
|----------|---|
| YKCA | Clinic Bundle A (1 each - YFBB, YFAA, YFAK, YFBI) |
| YKCB | Clinic Bundle B (1 each - YFAA, YFAK, YFBC, YFBE, YFBI, YFBJ, YFBL, YFBD) |
| YKCT | CT Positioning Sponge Bundle (1 each - YFBB, YFBC, YFBE, YFBF, YFBG, YFBI, YFBJ, YFBW) |
| YKGA | General Sponge Bundle A (1 each - YFAK, YFBC, YFBE, YFBI, YFBL) |
| YKGB | General Sponge Bundle B (1 each - YFAA, YFAK, YFBC, YFBE, YFBI, YFBJ, YFBL) |
| YKGC | General Sponge Bundle C (1 each - YFAA, YFAM, YFAR, YFAY, YFBL) |
| YKGD | General Sponge Bundle D (1 each - YFAL, YFAN, YFAR, YFBF, YFBJ, YFBL) |
| YKGG | General Sponge Bundle G (1 each - YFAA, YFAK, YFBD, YFBE, YFBI, YFBL) |
| YKGP | General Sponge Bundle P (1 each - YFAA, YFAK, YFBD, YFBI) |
| YKOA | Ortho Sponge Bundle A (1 each - YFAA, YFAM, YFAR, YFBD, YFBL) |
| YKPA | Pediatric Sponge Bundle (1 each - YFAC, YFAE, YFAK, YFAN, YFAQ, YFAX, YFBT) |
| YKWA | Wedge Sponge Bundle A (1 each - YFBC, YFBF, YFBG, YFBI) |
| YKWB | Wedge Sponge Bundle B (1 each - YFBD, YFBE, YFBF, YFBG) |
| YB-F-A02 | Urgent Care Sponge Bundle (1 each - YFAA, YFAD, YFAK, YFBC, YFBD, YFBE, YFBJ, YFBL, YFBN) |

VINYL SPONGES

**Measurements in L x W x H

RECTANGLE SPONGES

| |
|---|
| YVBR - Rectangle Sponge [6" x 4" x 1"] |
| YVBT - Rectangle Sponge [6" x 4" x 2"] |
| YVFG - Rectangle Bolster Sponge [12" x 8" x 6"] |
| YVFH - Rectangle Bolster Sponge [14" x 12" x 3"] |
| YVFJ - Rectangle Bolster Sponge [12" x 8" x 3"] |
| YVFK - Rectangle Bolster Sponge [20" x 13" x 3"] |
| YVFP - Rectangle Bolster Sponge [24" x 18" x 4"] |



YVFJ

WEDGE SPONGES

| | |
|--|--|
| YVBI - 15° Wedge Sponge Set [10" x 7" x 2"] | |
| YVDA - 45° Cardiac Wedge Sponge [28" x 14" x 14"] | |
| YVFO - Knee Wedge Bolster Sponge [17" x 13" x 6"] | YVFO |
| YVFC - 20° Wedge Bolster Sponge [24" x 18" x 6"] | YVFO - 45° Spinal Wedge Bolster Sponge [21" x 7" x 7"] |
| YVFD - 45° Wedge Bolster Sponge [16" x 8" x 8"] | YVFU - 30°- 60°- 90° Wedge Sponge [24" x 10" x 6"] |
| YVFE - 20° Incline Wedge Bolster Sponge [32" x 20" x 12"] | YVJV - Comfort Pillow Body Wedge [25" x 16" x 8.25"] |
| YVFF - 50° Incline Wedge Bolster Sponge [32" x 20" x 12"] | YVJZ - Endo-Ultrasound Wedge Sponge [20" x 20" x 6"] |
| YVFM - Knee Wedge Bolster Sponge [25" x 16" x 8.25"] | YVWA - Specialty Stretcher & Wheelchair Wedge Sponge [24.5" x 14" x 5.5-14"] |


SPECIALTY SPONGES



YVFA - Cylinder Bolster Sponge
[24" x 8" x 8"]



YVFB - Semi-Cylinder Bolster
Sponge [18" x 8" x 6"]



YVWS - Specialty Wheelchair Bolster
Sponge [17" x 13" x 6"]



YVAM - Flat Bottom V-Trough
[15" x 14" x 6.75"]



YVQ - DEXA Leg Bolster Sponge
[20" x 20" x 12"]



YB-V-A01 - Veterinary Sponge
Buncle

VINYL SANDBAGS & SANDSOCKS

Techno-Aide's premium sandbags and sandsocks are the ideal option for quick and comfortable patient positioning. These flexible positioners help position a patient's anatomy during difficult and specialty procedures. Easily make adjustments to reduce patient exposure and increase productivity. Techno-Aide's extensive line of positioning aids helps you save time, save costs and save lives!

SANDBAG BUNDLES & PAIRS

| | |
|-----------------|---|
| SBK-01-BLK-ST | Heavy Gauge Vinyl Sandbag General Bundle 1 - Standard Handles (Black) [3lb R (2x), 5lb R (1x), 5lb Sq (1x), 7lb R (2x), 10lb Sq (2x)] |
| SBK-01-BLU-ST | Heavy Gauge Vinyl Sandbag General Bundle 1 - Standard Handles (Blue) [3lb R (2x), 5lb R (1x), 5lb Sq (1x), 7lb R (2x), 10lb Sq (2x)] |
| SBK-02-BLK-ST | Heavy Gauge Vinyl Sandbag General Bundle 2 - Standard Handles (Black) [3lb R (2x), 7lb R (2x), 10lb R (2x)] |
| SBK-02-BLU-ST | Heavy Gauge Vinyl Sandbag General Bundle 2 - Standard Handles (Blue) [3lb R (2x), 7lb R (2x), 10lb R (2x)] |
| SBK-03-BLK-ST | Heavy Gauge Vinyl Sandbag Pediatric Bundle - Standard Handles (Black) [1lb R (1x), 2lb R (1x), 3lb R (1x), 5lb R (1x), 7lb R (1x)] |
| SBK-03-BLU-ST | Heavy Gauge Vinyl Sandbag Pediatric Bundle - Standard Handles (Blue) [1lb R (1x), 2lb R (1x), 3lb R (1x), 5lb R (1x), 7lb R (1x)] |
| SBK-04-BLK-ST | Heavy Gauge Vinyl Sandbag Femoral/Angio Bundle - Standard Handles (Black) [5lb R (1x), 7lb R (1x), 10lb R (1x)] |
| SBK-04-BLU-ST | Heavy Gauge Vinyl Sandbag Femoral/Angio Bundle - Standard Handles (Blue) [5lb R (1x), 7lb R (1x), 10lb R (1x)] |
| SBPR-05S-BLK-AC | Heavy Gauge Vinyl Sandbag Pair: Two 5lb (9x9) - AC Joint Handles (Black) |
| SBPR-05S-BLU-AC | Heavy Gauge Vinyl Sandbag Pair: Two 5lb (9x9) - AC Joint Handles (Blue) |
| SBPR-05S-BLK-ST | Heavy Gauge Vinyl Sandbag Pair: Two 5lb (9x9) - Standard Handles (Black) |
| SBPR-05S-BLU-ST | Heavy Gauge Vinyl Sandbag Pair: Two 5lb (9x9) - Standard Handles (Blue) |
| SBPR-10S-BLK-AC | Heavy Gauge Vinyl Sandbag Pair: Two 10lb (11x11) - AC Joint Handles (Black) |
| SBPR-10S-BLU-AC | Heavy Gauge Vinyl Sandbag Pair: Two 10lb (11x11) - AC Joint Handles (Blue) |
| SBPR-10S-BLK-ST | Heavy Gauge Vinyl Sandbag Pair: Two 10lb (11x11) - Standard Handles (Black) |
| SBPR-10S-BLU-ST | Heavy Gauge Vinyl Sandbag Pair: Two 10lb (11x11) - Standard Handles (Blue) |

VINYL SANDBAGS & SANDSOCKS


HEAVY GAUGE VINYL RECTANGULAR SANDBAGS: AC JOINT HANDLE

| | |
|--|--|
| SB-01R-BLK-AC - 1lb [4 x 5] SB-01R-BLU-AC - 1lb [4 x 5] |  <p>AC Joint Handle</p> |
| SB-02R-BLK-AC - 2lb [4 x 9] SB-02R-BLU-AC - 2lb [4 x 9] | |
| SB-03R-BLK-AC - 3lb [5 x 9] SB-03R-BLU-AC - 3lb [5 x 9] | |
| SB-05R-BLK-AC - 5lb [6 x 12] SB-05R-BLU-AC - 5lb [6 x 12] | |
| SB-07R-BLK-AC - 7lb [6 x 14] SB-07R-BLU-AC - 7lb [6 x 14] | |
| SB-10R-BLK-AC - 10lb [7 x 16] SB-10R-BLU-AC - 10lb [7 x 16] | |
| SB-12R-BLK-AC - 12lb [8 x 18] SB-12R-BLU-AC - 12lb [8 x 18] | |

HEAVY GAUGE VINYL RECTANGULAR SANDBAGS: STANDARD HANDLE

| | |
|--|--|
| SB-01R-BLK-ST - 1lb [4 x 5] SB-01R-BLU-ST - 1lb [4 x 5] |  <p>Standard Handle</p> |
| SB-02R-BLK-ST - 2lb [4 x 9] SB-02R-BLU-ST - 2lb [4 x 9] | |
| SB-03R-BLK-ST - 3lb [5 x 9] SB-03R-BLU-ST - 3lb [5 x 9] | |
| SB-05R-BLK-ST - 5lb [6 x 12] SB-05R-BLU-ST - 5lb [6 x 12] | |
| SB-07R-BLK-ST - 7lb [6 x 14] SB-07R-BLU-ST - 7lb [6 x 14] | |
| SB-10R-BLK-ST - 10lb [7 x 16] SB-10R-BLU-ST - 10lb [7 x 16] | |
| SB-12R-BLK-ST - 12lb [8 x 18] SB-12R-BLU-ST - 12lb [8 x 18] | |

HEAVY GAUGE VINYL SQUARE SANDBAGS: AC JOINT HANDLE

| | |
|--|--|
| SB-05S-BLK-AC - 5lb [9 x 9] SB-05S-BLU-AC - 5lb [9 x 9] |  <p>AC Joint Handle</p> |
| SB-07S-BLK-AC - 7lb [10 x 10] SB-07S-BLU-AC - 7lb [10 x 10] | |
| SB-10S-BLK-AC - 10lb [11 x 11] SB-10S-BLU-AC - 10lb [11 x 11] | |

HEAVY GAUGE VINYL SQUARE SANDBAGS: STANDARD HANDLE

| | |
|--|--|
| SB-05S-BLK-ST - 5lb [9 x 9] SB-05S-BLU-ST - 5lb [9 x 9] |  <p>Standard Handle</p> |
| SB-07S-BLK-ST - 7lb [10 x 10] SB-07S-BLU-ST - 7lb [10 x 10] | |
| SB-10S-BLK-ST - 10lb [11 x 11] SB-10S-BLU-ST - 10lb [11 x 11] | |

FLEXIBLE VINYL SANDSOCKS

| | |
|--|--|
| SS-01R-TAN-NH - 1lb [2 x 8] SS-01R-BLK-NH - 1lb [2 x 8] |  <p>Vinyl Sandsocks</p> |
| SS-03R-TAN-NH - 3lb [4 x 12] SS-03R-BLK-NH - 3lb [4 x 12] | |
| SS-04R-TAN-NH - 4lb [8 x 10] SS-04R-BLK-NH - 4lb [8 x 10] | SS-10R-TAN-NH - 10lb [6 x 20] SS-10R-BLK-NH - 10lb [6 x 20] |
| SS-05R-TAN-NH - 5lb [6 x 15] SS-05R-BLK-NH - 5lb [6 x 15] | SS-11R-TAN-NH - 11lb [6 x 28] SS-11R-BLK-NH - 11lb [6 x 28] |
| SS-07R-TAN-NH - 7lb [6 x 18] SS-07R-BLK-NH - 7lb [6 x 18] | SS-12R-TAN-NH - 12lb [6 x 25] SS-12R-BLK-NH - 12lb [6 x 25] |
| SS-08R-TAN-NH - 8lb [7 x 22] SS-08R-BLK-NH - 8lb [7 x 22] | SS-15R-TAN-NH - 15lb [6 x 40] SS-15R-BLK-NH - 15lb [6 x 40] |

FLEXIBLE VINYL SANDSOCKS GENERAL BUNDLES

| | |
|--|--|
| SSK-01-TAN-NH - Bundle 1 SSK-01-BLK-NH - Bundle 1 |  <p>SSK-01-TAN-NH</p> |
| SSK-02-TAN-NH - Bundle 2 SSK-02-BLK-NH - Bundle 2 | |
| SSK-03-TAN-NH - Bundle 3 SSK-03-BLK-NH - Bundle 3 | |

PLATFORM STOOLS & MRI POSITIONING AIDS

Techno-Aide offers an extensive line of positioning stools and specialty positioning supplies. Each item is carefully chosen to help get the right image the first time! Reduce patient exposure and increase productivity with these convenient positioning supplies and accessories. Techno-Aide's extensive line of positioning aids helps you save time, save costs, and save lives!



WMP-28



BBD-63

NYLON COVERED MRI POSITIONERS

MFP-01 - Rectangle Positioner Set (set of 2) [14" x 3" x 1"]

MFP-02 - Rectangle Positioner Set (set of 2) [5" x 3" x 1"]

MFP-03 - Wedge Positioner Set (set of 2) [10" x 10" x 3"]

MFP-04 - Wedge Positioner Set (set of 2) [14" x 4" x 1.5"]

MFP-05 - Rectangle Positioner Set (set of 2) [4" x 3" x 2"]

MFP-06 - Curved Wedge Positioner Set (set of 2) [5.2" x 5.3" x 3"]



MFP-01

MFP-07 - Semi-Round Neck Roll Positioner Set (set of 2) [7" x 3" x 2.3"]

MFP-10 - Knee Wedge Positioner Set (set of 2) [28.2" x 16" x 8"]

MFP-GNK - MRI Positioner Bundle (1 each - MFP-01, MFP-02, MFP-03, MFP-04, MFP-05, MFP-06)

MRI HEAVY-GAUGE VINYL SANDBAGS

MSB-NF03 - Square Sandbag: 3lb [8 x 8]

MSB-NF05 - Square Sandbag: 5lb [5 x 15]

MSB-NF07 - Square Sandbag: 7lb [9 x 12]

MSB-NF10 - Square Sandbag: 10lb [9 x 14]

MSB-NF12 - Square Sandbag: 12lb [12 x 12]

MSB-NFST - Sandbag Bundle - Standard Handle (1 each - MSB-NF03-NF12)



MSB-NF03



MSB-NFST - Sandbag Bundle

STANDARD & MRI PLATFORM STOOLS

| | |
|----------|--|
| BSH-65 | Standard Single Platform Step Stool with Handrail |
| BBS-61 | Single Platform Bariatric Step Stool with Handrail |
| BBD-63 | Dual Platform Bariatric Step Stool with Handrail |
| MBS-65NF | MRI Single Platform Step Stool with Handrail |
| MBD-68 | MRI Dual Platform Step Stool with Handrail |
| WMP-16 | Two Step Weight Bearing Imaging Platform (Comfort Closed Cell Steps) |
| WMP-22 | Two Step Weight Bearing Imaging Platform (Reinforced Plastic Steps) |
| WMP-14 | Three Step Weight Bearing Imaging Platform (Comfort Closed Cell Steps) |
| WMP-28 | Three Step Weight Bearing Imaging Platform (Reinforced Plastic Steps) |

SPECIALTY POSITIONING AIDS

Techno-Aide is proud to offer our extensive line of specialty positioning aids. This specially curated selection is designed to meet the needs of virtually any imaging facility, including MRI, CT, Veterinary, Fluoroscopy, and many others. Techno-Aide's specialty positioning aids are designed to save time, save costs, and save lives.



Patient Foot Immobilizer
FIL-12



Patient Head Immobilizer
RHI-AB

SPECIALTY POSITIONING AIDS

| | |
|------------|---|
| FIL-12 | Patient Foot Immobilizer |
| LGS | Adjustable Leg Support - Bariatric base [15" x 18"] (Previously LGS-B11) |
| LGS-S | Adjustable leg Support - Standard Base [9" x 12"] (Previously LGS-09) |
| RDS-12 | Universal Patient Restraint Strap |
| RHI-AB | Patient Head Immobilizer (table-Top Design with Coated Sponge Cushions) [Coated Sponge - 6 x 6 x 2] |
| RIS-60-BLK | Heavy-Gauge Vinyl Patient Immobilizer Strap [60" x 10"] |
| RIS-72-BLK | Heavy Gauge Vinyl Patient Immobilizer Strap [72" x 10"] |
| RIS-84-BLK | Heavy Gauge Vinyl Patient Immobilizer Strap [84" x 10"] |
| RIS-96-BLK | Heavy Gauge Vinyl Patient Immobilizer Strap [96" x 10"] |
| RLC-36 | Lateral Chest Image Support Positioner [36"] |
| RLC-48 | Lateral Chest Image Support Positioner [48"] |

Visit **TECHNO-AIDE.COM** to easily shop *all of our products* in one stop!



THE **ONE-STOP-SHOP**
FOR IMAGING SOLUTIONS

VETERINARY SPECIFIC POSITIONING AIDS

VETERINARY POSITIONING AIDS & BUNDLES


| | |
|------------|--|
| Q15645 | Thin Elastic Extremity Straps - 0.625 (5/8)" x 29" [Set of 2] |
| Q15647 | Standard Elastic Extremity Straps - 1" x 29" [Set of 2] |
| Q15648 | Extremity Straps Table Liner- 15' (5 yds) [Black] |
| Q15645-B | Thin Elastic Extremity Straps Bundle - 0.625 (5/8)" x 29" [Set of 2] + Adhesive Table Liner [15' Black Loop] |
| Q15647-B | Standard Elastic Extremity Straps Bundle - 1" x 29" [Set of 2] + Adhesive Table Liner [15' Black Loop] |
| VIS-60-BLK | Heavy-Gauge Vinyl Veterinary Immobilizer Strap (60" x 5") |
| VIS-72-BLK | Heavy-Gauge Vinyl Veterinary Immobilizer Strap (72" x 5") |
| VIS-84-BLK | Heavy-Gauge Vinyl Veterinary Immobilizer Strap (84" x 5") |
| VIS-96-BLK | Heavy-Gauge Vinyl Veterinary Immobilizer Strap (96" x 5") |
| VIT-B | Heavy-Gauge Veterinary Vinyl Immobilizer Trough - Bundle (one of each size) |
| VIT-L | Heavy-Gauge Veterinary Vinyl Immobilizer Trough - Large (24" x 11" x 7") |
| VIT-M | Heavy-Gauge Veterinary Vinyl Immobilizer Trough - Medium (15" x 7" x 5") |
| VIT-S | Heavy-Gauge Veterinary Vinyl Immobilizer Trough - Small (10" x 4" x 3") |
| VIT-X | Heavy-Gauge Veterinary Vinyl Immobilizer Trough - XL (36" x 14" x 9") |
| YCCA | Coated Small Extremity Sponge [12.5" x 7" x 3"] |
| YCCB | Coated Medium Extremity Sponge [19.75" x 11.5" x 5"] |
| YCCC | Coated Large Extremity Sponge [30.5" x 11" x 7"] |
| YCCD | Coated XL Extremity Sponge [44.5" x 14" x 9"] |
| YFCA | Non-Coated Coated Small Extremity Sponge [12.5" x 7" x 3"] |
| YFCB | Non-Coated Medium Extremity Sponge [19.75" x 11.5" x 5"] |
| YFCC | Non-Coated Large Extremity Sponge [24" x 11" x 7"] |
| YFCD | Non-Coated XL Extremity Sponge [36" x 14" x 9"] |

COMING SOON!

| | |
|------|--|
| YXCA | Closed Cell Small Extremity Sponge [13" x 7" x 3.5"] |
| YXCB | Closed Cell Medium Extremity Sponge [19.5" x 11.5" x 5"] |
| YXCC | Closed Cell Large Extremity Sponge [30.5" x 11" x 7"] |
| YXCD | Closed Cell XL Extremity Sponge [44" x 14" x 9"] |

STANDARD TABLE PADS

- Comfortable and durable vinyl over a premium foam interior
- Perfect for x-ray tables and stretchers
- Increases patient comfort and stability
- Available in universal sizes designed to fit virtually every standard make and model table
- Feature tapered ends with brass grommets for proper hanging storage
- Made from radiolucent foam and material
- Easy to clean and sanitize!
- Available in three premium foam types
- Choose from black or blue high-quality vinyl, or our deluxe reinforced Techno-Tuff™ vinyl



Standard Radiolucent X-Ray Table Pads

PREMIUM TABLE PADS

- Upgraded foam interior options, including a deluxe memory foam style core
- Perfect for x-ray tables and stretchers
- Increases patient comfort and stability
- Available in universal sizes designed to fit virtually every standard make and model table
- Feature sleek sewn and sealed boxed corners
- Made from radiolucent foam and material
- Easy to clean and sanitize!
- Available in three premium foam types



Premium Radiolucent X-Ray Table Pads

ECOMONY TABLE PADS

- Perfect for x-ray tables and stretchers
- Increases patient comfort and stability
- Available in universal sizes designed to fit virtually every standard make and model table
- Feature tapered ends with brass grommets for proper hanging storage
- Made from radiolucent foam and material
- Easy to clean and sanitize!
- Available in three premium foam types



Premium Radiolucent X-Ray Table Pads

WALL MOUNTED TABLE PAD STORAGE RACK

Designed to work specifically with Techno-Aide grommated table pads, this durable low-profile wall mounted rack helps extend the life of your table pads.



PAD-R2